

## LogixProgram L27 Ver4.ACD

pdo

Controller test

Controller Fault Handler

Power-Up Handler

## Tasks

## (P) P200ms\_p10

GUI

Handle\_GUI

ValveProfile

Main

## (P) P\_Sample\_p5

Get raw analog input and control of the breathing - MUST HAVE HIGHEST priority of all peridoc tasks

measArray

Main

FirstScan

measArrays

ISR\_ValveCtrl

Main

## (P) P\_Sample\_p8

Setup

Main

InitGUI

InitISR\_RS232

InitVariables

calcBreath

Main

initVent

Main

initISR\_ValveCtrl

Init\_ISR\_PressWatchID

BreathDetect

Main

measCalc

Main

MainFlow

Main

Unscheduled

## Motion Groups

Ungrouped Axes

DymmyMaster

DymmySlave

## Add-On Instructions

AOI\_ARRAYSHIFT

Logic

AOI\_Ceil

Logic

AOI\_FillArrayInt

Logic

AOI\_FillArrayReal

Logic

AOI\_FillArrayReal50

Logic

AOI\_FillArray\_Real2D

Logic

AOI\_FilterSignal

Logic

Prescan

AOI\_Floor

Logic

AOI\_Max

Logic

- Prescan
- AOI\_Mean2D
  - Logic
  - Prescan
- AOI\_Mean50
  - Logic
  - Prescan
- AOI\_Round
  - Logic
- AOI\_Trunc
  - Logic

## Data Types

User-Defined

Strings

Add-On-Defined

- AOI\_ARRAYSHIFT
- AOI\_Ceil
- AOI\_FillArrayInt
- AOI\_FillArrayReal
- AOI\_FillArrayReal50
- AOI\_FillArray\_Real2D
- AOI\_FilterSignal
- AOI\_Floor
- AOI\_Max
- AOI\_Mean2D
- AOI\_Mean50
- AOI\_Round
- AOI\_Trunc

Module-Defined

- <sup>101</sup><sub>010</sub> AB:1769\_HSC1\_Range:C:0
- <sup>101</sup><sub>010</sub> AB:1769\_IF4FXOF2F:C:0
- <sup>101</sup><sub>010</sub> AB:1769\_IF4FXOF2F:I:0
- <sup>101</sup><sub>010</sub> AB:1769\_IF4FXOF2F:O:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_AnalogIO1:C:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_AnalogIO1:I:0<sup>101</sup><sub>010</sub>
- AB:Embedded\_AnalogIO1:O:0<sup>101</sup><sub>010</sub>
- AB:Embedded\_DiscreteIO1:C:0<sup>101</sup><sub>010</sub>
- AB:Embedded\_DiscreteIO1:I:0<sup>101</sup><sub>010</sub>
- AB:Embedded\_DiscreteIO1:O:0<sup>101</sup><sub>010</sub>
- AB:Embedded\_DiscreteIO:C:0<sup>101</sup><sub>010</sub>
- AB:Embedded\_DiscreteIO:I:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_DiscreteIO:O:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_HSC1:C:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_HSC1:I:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_HSC1:O:0
- <sup>101</sup><sub>010</sub> AB:Embedded\_HSC1\_STRUCT\_OUT1:O:0

## Trends

CalcBreath  
 Test

## I/O Configuration

1769 Bus

- [0] 1769-L27ERM-QBFC1B test
  - [1] Embedded Discrete\_IO
  - [2] Embedded Analog\_IO
  - [3] Embedded Counters

Ethernet

- 1769-L27ERM-QBFC1B test

**General**

Vendor:	Rockwell Automation/Allen-Bradley	Mode:	Offline
Revision:	32.013	Key Switch Position:	Offline
Chassis Type:	<none>	Created:	26-03-2020 15:02:52
Slot:	0	Edited:	02-04-2020 20:43:03

**Date/Time**

Date and Time:	<offline>
Time Zone:	<offline>
Daylight Saving (+00:00):	<offline>
Enable Time Synchronization:	No
Is the system time master:	<offline>
Is a synchronized time slave:	<offline>
Duplicate CST Master Detected:	<offline>
CST Mastership disabled:	<offline>
No CST Master:	<offline>

**Advanced Time Sync**

CIP Sync Time Synchronization:	Disabled
--------------------------------	----------

**Advanced**

Controller Fault Handler:	<none>	Match Project To Controller:	No
Power-Up Handler:	<none>	Serial Number:	C01FB861
System Overhead Time Slice:	20 %	Allow Consumed Tags to Use RPI Provided by Producer:	No

During unused System Overhead Time

Slice: Run Continuous Task

**SFC Execution**

Execution Control:	Execute current active steps only	Last Scan of Active Step:	Don't scan
Restart Position:	Restart at most recently executed step		

**Nonvolatile Memory**

&lt;offline&gt;

**Memory (Estimate)**

Memory Option:	1769-L27ERM-QBFC1B		
Estimated I/O Memory			
Total Memory:	1,048,576 bytes	Max Used:	4,640 bytes
Free Memory:	1,043,936 bytes	Largest Block Free:	1,043,936 bytes
Used Memory:	4,640 bytes		
Estimated Data and Logic Memory			
Total Memory:	1,048,576 bytes	Max Used:	665,744 bytes
Free Memory:	382,832 bytes	Largest Block Free:	382,832 bytes
Used Memory:	665,744 bytes		

**Security:**

Primary Security Authority:	No Protection	Restrict Communications Except Through Selected Slots:	No
Use only the selected Security Authority for Authentication and Authorization:	No		
Secondary Security Authority:	No Protection	Changes To Detect:	16#ffff_ffff_ffff_ffff
Use only the selected Security Authority for Authentication and Authorization:	No	Audit Value:	<offline>
Permission Set:	No		

**Port Configuration****Port 1**

Enable:

Yes

**Port 2**

Enable:

Yes

## LogixProgram\_L27\_Ver4.ACD

Name	Value	Data Type	Scope
_CopDes1	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopDes1 - ISR_ValveCtrl/Main - #142, #143, #144, *#140			
_CopDes4	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopDes4 - ISR_ValveCtrl/Main - #163, #164, #165, #168, *#161			
_CopDes5	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopDes5 - ISR_ValveCtrl/Main - #170, #171, #172, *#168			
_CopDes6	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopDes6 - ISR_ValveCtrl/Main - #179, #180, #181, *#177			
_CopDes8	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopDes8 - ISR_ValveCtrl/Main - #198, #199, #200, *#196			
_CopLength	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength - ISR_ValveCtrl/Main - #135, #136, #137, #140, *#134			
_CopLength1	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength1 - ISR_ValveCtrl/Main - #142, #143, #144, *#141			
_CopLength3	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength3 - ISR_ValveCtrl/Main - #150, #151, #152, *#149			
_CopLength4	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength4 - ISR_ValveCtrl/Main - #163, #164, #165, *#162			
_CopLength5	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength5 - ISR_ValveCtrl/Main - #170, #171, #172, *#169			
_CopLength6	0	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength6 - ISR_ValveCtrl/Main - #179, #180, #181, *#178			
_CopLength7	2060	DINT	test
Constant	No		
External Access:	Read/Write		
_CopLength7 - ISR_ValveCtrl/Main - #190, #191, #192, *#189			
_CopLength8	0	DINT	test
Constant	No		
External Access:	Read/Write		

## LogixProgram\_L27\_Ver4.ACD

<b>_CopLength8 (Continued)</b>			
<i>_CopLength8 - ISR_ValveCtrl/Main</i>	- #198, #199, #200, *#197		
<b>_CopSource</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource - ISR_ValveCtrl/Main</i>	- #135, #136, #137, *#133		
<b>_CopSource1</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource1 - ISR_ValveCtrl/Main</i>	- #142, #143, #144, *#139		
<b>_CopSource3</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource3 - ISR_ValveCtrl/Main</i>	- #150, #151, #152, *#148		
<b>_CopSource4</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource4 - ISR_ValveCtrl/Main</i>	- #163, #164, #165, *#160		
<b>_CopSource5</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource5 - ISR_ValveCtrl/Main</i>	- #170, #171, #172, *#167		
<b>_CopSource6</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource6 - ISR_ValveCtrl/Main</i>	- #179, #180, #181, *#176		
<b>_CopSource7</b>	-30	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource7 - ISR_ValveCtrl/Main</i>	- #190, #191, #192, *#188		
<b>_CopSource8</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>_CopSource8 - ISR_ValveCtrl/Main</i>	- #198, #199, #200, *#195		
<b>_edbAoiFilterFlowAir</b>		AOI_FilterSignal	test
Constant	No		
External Access:	Read/Write		
<i>_edbAoiFilterFlowAir - measArray/measArrays</i>	- *#109		
<b>_edbAoiFilterFlowAir.EnableIn</b>	1	BOOL	
Enable Input - System Defined Parameter			
<b>_edbAoiFilterFlowAir.EnableOut</b>	1	BOOL	
Enable Output - System Defined Parameter			
<b>_edbAoiFilterFlowAir.In_bInitialize</b>	0	BOOL	
<i>_edbAoiFilterFlowAir.In_bInitialize - measArray/measArrays</i>	- *#101		
<b>_edbAoiFilterFlowAir.Par_eType</b>	0	DINT	
0: MEDIAN; 1: MEAN			
<i>_edbAoiFilterFlowAir.Par_eType - measArray/measArrays</i>	- *#84		
<b>_edbAoiFilterFlowAir.Par_nWindowLeft</b>	5	DINT	
<i>_edbAoiFilterFlowAir.Par_nWindowLeft - measArray/measArrays</i>	- *#82		
<b>_edbAoiFilterFlowAir.Par_nWindowsRight</b>	5	DINT	
<i>_edbAoiFilterFlowAir.Par_nWindowsRight - measArray/measArrays</i>	- *#83		
<b>_edbAoiFilterFlowAir.Par_bUseWeights</b>			

## LogixProgram\_L27\_Ver4.ACD

<b>_edbAoiFilterFlowAir (Continued)</b>	0	BOOL	
Used for mean filtering			
<i>_edbAoiFilterFlowAir.Par_bUseWeights - measArray/measArrays - *#85</i>			
<b>_edbAoiFilterFlowAir.Out_Buffering</b>	0	BOOL	
<i>_edbAoiFilterFlowAir.Out_Buffering - measArray/measArrays - #116</i>			
<b>_edbAoiFilterFlowAir.Out_rFilteredSignalDelayed</b>	980.0	REAL	
<i>_edbAoiFilterFlowAir.Out_rFilteredSignalDelayed - measArray/measArrays - #122</i>			
<b>_edbAoiFilterFlowO2</b>		AOI_FilterSignal	test
Constant	No		
External Access:	Read/Write		
<i>_edbAoiFilterFlowO2 - measArray/measArrays - *#111</i>			
<b>_edbAoiFilterFlowO2.EnableIn</b>	1	BOOL	
Enable Input - System Defined Parameter			
<b>_edbAoiFilterFlowO2.EnableOut</b>	1	BOOL	
Enable Output - System Defined Parameter			
<b>_edbAoiFilterFlowO2.In_bInitialize</b>	0	BOOL	
<i>_edbAoiFilterFlowO2.In_bInitialize - measArray/measArrays - *#102</i>			
<b>_edbAoiFilterFlowO2.Par_eType</b>	0	DINT	
0: MEDIAN; 1: MEAN			
<i>_edbAoiFilterFlowO2.Par_eType - measArray/measArrays - *#89</i>			
<b>_edbAoiFilterFlowO2.Par_nWindowLeft</b>	5	DINT	
<i>_edbAoiFilterFlowO2.Par_nWindowLeft - measArray/measArrays - *#87</i>			
<b>_edbAoiFilterFlowO2.Par_nWindowsRight</b>	5	DINT	
<i>_edbAoiFilterFlowO2.Par_nWindowsRight - measArray/measArrays - *#88</i>			
<b>_edbAoiFilterFlowO2.Par_bUseWeights</b>	0	BOOL	
Used for mean filtering			
<i>_edbAoiFilterFlowO2.Par_bUseWeights - measArray/measArrays - *#90</i>			
<b>_edbAoiFilterFlowO2.Out_Buffering</b>	0	BOOL	
<i>_edbAoiFilterFlowO2.Out_Buffering - measArray/measArrays - #117</i>			
<b>_edbAoiFilterFlowO2.Out_rFilteredSignalDelayed</b>	981.0	REAL	
<i>_edbAoiFilterFlowO2.Out_rFilteredSignalDelayed - measArray/measArrays - #126</i>			
<b>_edbAoiFilterPress</b>		AOI_FilterSignal	test
Constant	No		
External Access:	Read/Write		
<i>_edbAoiFilterPress - measArray/measArrays - *#110</i>			
<b>_edbAoiFilterPress.EnableIn</b>	1	BOOL	
Enable Input - System Defined Parameter			
<b>_edbAoiFilterPress.EnableOut</b>	1	BOOL	
Enable Output - System Defined Parameter			
<b>_edbAoiFilterPress.In_bInitialize</b>	0	BOOL	
<i>_edbAoiFilterPress.In_bInitialize - measArray/measArrays - *#103</i>			
<b>_edbAoiFilterPress.Par_eType</b>	0	DINT	
0: MEDIAN; 1: MEAN			
<i>_edbAoiFilterPress.Par_eType - measArray/measArrays - *#94</i>			
<b>_edbAoiFilterPress.Par_nWindowLeft</b>	5	DINT	
<i>_edbAoiFilterPress.Par_nWindowLeft - measArray/measArrays - *#92</i>			
<b>_edbAoiFilterPress.Par_nWindowsRight</b>	5	DINT	
<i>_edbAoiFilterPress.Par_nWindowsRight - measArray/measArrays - *#93</i>			
<b>_edbAoiFilterPress.Par_bUseWeights</b>	0	BOOL	
Used for mean filtering			

## LogixProgram\_L27\_Ver4.ACD

<b>_edbAoiFilterPress (Continued)</b>			
<code>_edbAoiFilterPress.Par_bUseWeights - measArray/measArrays - *#95</code>			
<b>_edbAoiFilterPress.Out_Buffering 0</b>	BOOL		
<code>_edbAoiFilterPress.Out_Buffering - measArray/measArrays - #118</code>			
<b>_edbAoiFilterPress.Out_rFilteredSignalDelayed 3947.0</b>	REAL		
<code>_edbAoiFilterPress.Out_rFilteredSignalDelayed - measArray/measArrays - #124</code>			
<b>_edbCanStore 1</b>	BOOL		test
Constant No			
External Access: Read/Write			
<code>_edbCanStore - measArray/measArrays - #120, *#116</code>			
<b>_edbFlowArrFiltered</b>	REAL[1000,2]		test
Constant No			
External Access: Read/Write			
<code>_edbFlowArrFiltered - MainFlow/Main - *#62</code>			
<code>_edbFlowArrFiltered - measArray/measArrays - *#122, *#123</code>			
<b>_edbFlowO2Filtered</b>	REAL[1000,2]		test
Constant No			
External Access: Read/Write			
<code>_edbFlowO2Filtered - MainFlow/Main - *#63</code>			
<code>_edbFlowO2Filtered - measArray/measArrays - *#124, *#125</code>			
<b>_edbHalfWindowSize 5</b>	DINT		test
Constant No			
External Access: Read/Write			
<code>_edbHalfWindowSize - measArray/measArrays - #119, #121, #82, #83, #87, #88, #92, #93, *#67</code>			
<b>_edbIdx 590</b>	DINT		test
Constant No			
External Access: Read/Write			
<code>_edbIdx - measArray/measArrays - #122, #123, #124, #125, #126, #127, *#121</code>			
<b>_edbParUseWeights 0</b>	BOOL		test
Constant No			
External Access: Read/Write			
<code>_edbParUseWeights - measArray/measArrays - #85, #90, #95, *#74</code>			
<b>_edbParWeights</b>	REAL[31]		test
Constant No			
External Access: Read/Write			
<code>_edbParWeights - measArray/measArrays - *#109, *#110, *#111</code>			
<b>_edbPressArrFiltered</b>	REAL[1000,2]		test
Constant No			
External Access: Read/Write			
<code>_edbPressArrFiltered - MainFlow/Main - *#64</code>			
<code>_edbPressArrFiltered - measArray/measArrays - *#126, *#127</code>			
<b>_edbTypeOffFilter 0</b>	DINT		test
Constant No			
External Access: Read/Write			
<code>_edbTypeOffFilter - measArray/measArrays - #84, #89, #94, *#73</code>			
<b>_endPosition 0</b>	DINT		test
Constant No			
External Access: Read/Write			
<code>_endPosition - ISR_ValveCtrl/Main - #104, #106, #121, #125, #127, #156, #194, *#105, *#121</code>			
<code>_endPosition - measArray/FirstScan - *#107</code>			
<b>_endPositionOld 0</b>	DINT		test
Constant No			

## LogixProgram\_L27\_Ver4.ACD

**\_endPositionOld (Continued)**

External Access: Read/Write  
*\_endPositionOld - ISR\_ValveCtrl/Main - #104, \*#106, \*#205*

**\_EndSendPosition**

0 DINT test  
Constant No  
External Access: Read/Write  
*\_EndSendPosition - ISR\_ValveCtrl/Main - #158, #169, #178, #197, \*#156, \*#194*

**\_i**

1000 DINT test  
Constant No  
External Access: Read/Write  
*\_i - measCalc/Main - #12, #18, #18, #19, #27, #27, #29, #30, #32, #41, #41, #42, \*#10*

**\_J**

196 INT test  
Constant No  
External Access: Read/Write  
*\_J - measCalc/Main - #108, \*#107*

**\_n**

412 DINT test  
Constant No  
External Access: Read/Write  
*\_n - measCalc/Main - \*#139*

**\_PointAfterEnd**

30 DINT test  
Constant No  
External Access: Read/Write  
*\_PointAfterEnd - ISR\_ValveCtrl/Main - #122, #125, #156, #194, \*#122*  
*\_PointAfterEnd - measArray/FirstScan - \*#108*

**\_StartPointBerfore**

30 DINT test  
Constant No  
External Access: Read/Write  
*\_StartPointBerfore - ISR\_ValveCtrl/Main - #119, #125, #129, #149, #161, #177, #187, \*#119*  
*\_StartPointBerfore - measArray/FirstScan - \*#105*

**\_StartPosition**

0 DINT test  
Constant No  
External Access: Read/Write  
*\_StartPosition - ISR\_ValveCtrl/Main - #120, #125, #127, #129, #160, #162, #168, #176, #178, #187, #28, #30, \*#120, \*#29*  
*\_StartPosition - measArray/FirstScan - \*#106*

**\_StartpositionOld**

0 DINT test  
Constant No  
External Access: Read/Write  
*\_StartpositionOld - ISR\_ValveCtrl/Main - #28, \*#206, \*#30*

**\_T**

0 DINT test  
Constant No  
External Access: Read/Write  
*\_T - measCalc/Main - #149, #149, #150, #150, #151, #151, #152, #152, #153, #153, #154, #154, #155, #155, #156, #156, #157, #157, #158, #158, #159, #159, \*#147*

**\_z**

51 DINT test  
Constant No  
External Access: Read/Write  
*\_z - measCalc/Main - #183, #183, #187, #188, #194, #195, #198, \*#182, \*#185, \*#190, \*#199*

**AirConstantFlow**

68.35443 REAL test  
Constant No  
External Access: Read/Write  
*AirConstantFlow - calcBreath/Main - #46, #66, \*#41*

**AirDown**

CAM\_PROFILE[1010]

test

## LogixProgram\_L27\_Ver4.ACD

**AirDown (Continued)**

Constant No  
External Access: Read/Write  
*AirDown - ValveProfile/Main - 0(MAPC)*

**AirmADer** 0.0 REAL test

Constant No  
External Access: Read/Write  
*AirmADer - ISR\_ValveCtrl/Main - \*#63*

**AirmAOut** 0.0 REAL test

Constant No  
External Access: Read/Write  
*AirmAOut - ISR\_ValveCtrl/Main - #64, \*#63*

**AirmASlope** 49211.336 REAL test

Constant No  
External Access: Read/Write  
*AirmASlope - ISR\_ValveCtrl/Main - \*#63*

**AirRampFlowIncreasePerAdjust** 3.4177215 REAL test

Constant No  
External Access: Read/Write  
*AirRampFlowIncreasePerAdjust - calcBreath/Main - #63, \*#46*

**AirUp** CAM\_PROFILE[1010] test

Constant No  
External Access: Read/Write  
*AirUp - ISR\_ValveCtrl/Main - #63*  
*AirUp - ValveProfile/Main - 0(MAPC)*

**AnalogInputFlow** 1059 INT test

AliasFor: Local:2:I.Ch3Data  
Base Tag: Local:2:I.Ch3Data  
Constant No  
External Access: Read/Write  
*AnalogInputFlow - measArray/measArrays - #32*

**AnalogInputFlowAir** 980 INT test

AliasFor: Local:2:I.Ch1Data  
Base Tag: Local:2:I.Ch1Data  
Constant No  
External Access: Read/Write  
*AnalogInputFlowAir - measArray/measArrays - #109, #34*

**AnalogInputFlowO2** 981 INT test

AliasFor: Local:2:I.Ch0Data  
Base Tag: Local:2:I.Ch0Data  
Constant No  
External Access: Read/Write  
*AnalogInputFlowO2 - measArray/measArrays - #111, #36*

**AnalogInputPressure** 3943 INT test

AliasFor: Local:2:I.Ch2Data  
Base Tag: Local:2:I.Ch2Data  
Constant No  
External Access: Read/Write  
*AnalogInputPressure - ISR\_ValveCtrl/Main - #1*  
*AnalogInputPressure - measArray/measArrays - #110, #38*

**AnalogOutFlowAir** 0 INT test

AliasFor: Local:2:O.Ch1Data  
Base Tag: Local:2:O.Ch1Data  
Constant No

LogixProgram\_L27\_Ver4.ACD

**AnalogOutFlowAir (Continued)**

External Access: Read/Write  
*AnalogOutFlowAir - ISR\_ValveCtrl/Main - \*#53, \*#64, \*#99*

**AnalogOutFlowO2** 0 INT test

AliasFor: Local:2:O.Ch0Data  
Base Tag: Local:2:O.Ch0Data  
Constant No  
External Access: Read/Write  
*AnalogOutFlowO2 - ISR\_ValveCtrl/Main - \*#52, \*#62, \*#98*

**avgCrsDyn** 0.44807863 REAL test

1 min avg of dynamic respiratory system compliance in ml/cmH2O, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgCrsDyn - measCalc/Main - \*#221*

**avgCrsQStat** 0.50136733 REAL test

1 min avg of quasistatic respiratory system compliance in ml/cmH2O, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgCrsQStat - measCalc/Main - \*#223*

**avgFlowImax** 18987.0 REAL test

1 min avg of max insp flow in l/min, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgFlowImax - measCalc/Main - \*#211*

**avgMV** 41214.273 REAL test

1 min avg of minute ventilation in l/min, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgMV - measCalc/Main - \*#219*

**avgPEEP** 4402.0 REAL test

1 min avg of positive end expiratory pressure in cmH2O, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgPEEP - measCalc/Main - \*#209*

**avgPpeak** 10534.0 REAL test

1 min avg of peak insp pressure in cmH2O, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgPpeak - measCalc/Main - \*#207*

**avgPplat** 9882.25 REAL test

1 min avg of plateau insp pressure in cmH2O, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgPplat - measCalc/Main - \*#213*

**avgRaw** 2.059567 REAL test

1 min avg of airway resistance in cmH2O/l/s, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgRaw - measCalc/Main - \*#225*

**avgRR** 15.0 REAL test

1 min avg of respiratory rate in 1/min, from measAvgCalc()  
Constant No  
External Access: Read/Write  
*avgRR - measCalc/Main - \*#215*

LogixProgram\_L27\_Ver4.ACD

<b>avgVti</b>	2747.6182	REAL	test
	1 min avg of insp tidal volume in ml, from measAvgCalc()		
Constant	No		
External Access:	Read/Write		
<i>avgVti - measCalc/Main - *#217</i>			
<b>BackingTag</b>		MOTION_INSTRUCTION	test
Constant	No		
External Access:	Read/Write		
<i>BackingTag - ValveProfile/Main - *0(MAPC)</i>			
<b>BreathDuration</b>	4	INT	test
Constant	No		
External Access:	Read/Write		
<i>BreathDuration - measCalc/Main - #I73, #65, *#64</i>			
<b>BreathDurationArr</b>		REAL[50]	test
50 length array with seconds for each breath for calculation of minute avgVentStateArr. Calcualted from knowing meas at 100Hz for all breath data			
Constant	No		
External Access:	Read/Write		
<i>BreathDurationArr - initVent/initISR_ValveCtrl - *#5</i>			
<i>BreathDurationArr - measArray/FirstScan - *#52</i>			
<i>BreathDurationArr - measCalc/Main - #I59, #I83, #I88, #I95, *#I59</i>			
<i>BreathDurationArr[0] - measCalc/Main - *#173</i>			
<b>BreathDurationInMS</b>	4000.0002	REAL	test
Constant	No		
External Access:	Read/Write		
<i>BreathDurationInMS - calcBreath/Main - #11, #12, #15, #20, #21, *#10</i>			
<b>CalcNewData</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>CalcNewData - ISR_ValveCtrl/Main - *#207, *#211</i>			
<i>CalcNewData - measCalc/Main - #4, *#227</i>			
<b>ConstansFlowNrCalc</b>	40.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>ConstansFlowNrCalc - calcBreath/Main - #49, *#48</i>			
<b>ConstantFlowNrOfAdjustments</b>	40.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>ConstantFlowNrOfAdjustments - calcBreath/Main - #64, #67, #67, *#49</i>			
<b>ConstantFlowTimeInMS</b>	400.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>ConstantFlowTimeInMS - calcBreath/Main - #27, #48, *#22</i>			
<b>ConstantInspiratoryFlow</b>	72.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>ConstantInspiratoryFlow - calcBreath/Main - #31, #40, #41, *#27</i>			
<b>counter</b>	594	DINT	test
Counter for Flow fill array			
Constant	No		
External Access:	Read/Write		
<i>counter - ISR_ValveCtrl/Main - #105, #29</i>			
<i>counter - measArray/FirstScan - *#3</i>			

## LogixProgram\_L27\_Ver4.ACD

**counter (Continued)***counter - measArray/measArrays - #119, #121, #132, #135, #31, #31, #32, #33, #34, #35, #36, #37, #38, #39, #58, \*#132, \*#136*

<b>CrsDyn</b>	0.5221349	REAL	test
dynamic respiratory system compliance in ml/cmH2O			
Constant	No		
External Access:	Read/Write		

*CrsDyn - measCalc/Main - #170, \*#127*

<b>CrsDynArr</b>	REAL[50]	test
50 length array with dynamic respiratory system compliance in ml/cmH2O, from measCalc()		
Constant	No	
External Access:	Read/Write	
<i>CrsDynArr - measArray/FirstScan - *#49</i>		
<i>CrsDynArr - measCalc/Main - #156, *#156, *#220</i>		
<i>CrsDynArr[0] - measCalc/Main - *#170</i>		

<b>CrsQStat</b>	0.59738743	REAL	test
quasistatic respiratory system compliance in ml/cmH2O			
Constant	No		
External Access:	Read/Write		

*CrsQStat - measCalc/Main - #171, \*#131, \*#133*

<b>CrsQStatArr</b>	REAL[50]	test
50 length array with quasistatic respiratory system compliance in ml/cmH2O, from measCalc()		
Constant	No	
External Access:	Read/Write	
<i>CrsQStatArr - measArray/FirstScan - *#50</i>		
<i>CrsQStatArr - measCalc/Main - #157, *#157, *#222</i>		
<i>CrsQStatArr[0] - measCalc/Main - *#171</i>		

<b>CsumData</b>	659	DINT	test
Constant	No		
External Access:	Read/Write		
<i>CsumData - MainFlow/Main - #11, #12, *#9</i>			

<b>CsumDataOld</b>	659	DINT	test
Constant	No		
External Access:	Read/Write		
<i>CsumDataOld - MainFlow/Main - #11, *#12</i>			

<b>data_changed</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>data_changed - GUI/Handle_GUI - #139, #194, #31, #43, #92, *#140, *#143, *#185, *#44, *#48, *#93, *#96</i>			
<i>data_changed - measArray/FirstScan - *#28</i>			

<b>DymmyMaster</b>		AXIS_VIRTUAL	test
External Access:	Read/Write		

*DymmyMaster - ValveProfile/Main - 0(MAPC)*

<b>DymmySlave</b>		AXIS_VIRTUAL	test
External Access:	Read/Write		

*DymmySlave - ValveProfile/Main - 0(MAPC)*

<b>endI</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>endI - measArray/FirstScan - *#39</i>			
<i>endI - measCalc/Main - #206, #208, #210, #212, #214, #216, #218, #220, #222, #224, *#189, *#198</i>			

<b>EndPoint</b>	194	INT	test
Constant	No		
External Access:	Read/Write		

LogixProgram\_L27\_Ver4.ACD**EndPoint (Continued)***EndPoint - measCalc/Main - #99, \*#98*

Exp	0	BOOL	test
-----	---	------	------

Constant No  
External Access: Read/Write

*Exp - ISR\_ValveCtrl/Main - #34, #6, #6, \*#115, \*#13, \*#17, \*#72, \*#8*

*Exp - MainFlow/Main - #53, \*#67*

expCount	0	INT	test
----------	---	-----	------

between 1-1000 to indicate where valveCtrl should read from expOpenCloseCtrl[]

Constant No  
External Access: Read/Write

*expCount - ISR\_ValveCtrl/Main - #100, #102, #15, #7, #96, \*#42, \*#96*

*expCount - MainFlow/Main - #54, \*#68*

expCountLength	300	INT	test
----------------	-----	-----	------

between 1-1000 to indicate when valveCtrl should stop reading from expOpenCloseCtrl[], copied from nextExpCountLength by ISR\_ValveCtrl()

Constant No  
External Access: Read/Write

*expCountLength - initVent/Main - #11, \*#9*

*expCountLength - ISR\_ValveCtrl/Main - #102, #15, #7, \*#224*

*expCountLength - MainFlow/Main - #54*

expCountlengthCalc	300.00003	REAL	test
--------------------	-----------	------	------

Constant No  
External Access: Read/Write

*expCountlengthCalc - calcBreath/Main - #17, \*#16*

expCountLengthNew	300.0	REAL	test
-------------------	-------	------	------

Constant No  
External Access: Read/Write

*expCountLengthNew - calcBreath/Main - #74, #83, \*#17*

expEndI	0	INT	test
---------	---	-----	------

Constant No  
External Access: Read/Write

*expEndI - measArray/FirstScan - \*#37*

*expEndI - measCalc/Main - #139, #37, #57, #64, \*#42, \*#58, \*#8*

expOpen	1	BOOL	test
---------	---	------	------

AliasFor: Local:1:O.Data.0  
Base Tag: Local:1:O.Data.0  
Constant No  
External Access: Read/Write

*expOpen - ISR\_ValveCtrl/Main - \*#100, \*#54, \*#60*

expOpenCloseCtrl	INT[1000]	test
------------------	-----------	------

1000 length array with open (1) or close (0) for exp valve, moved from nextexpOpenCloseCtrl[] by ISR\_ValveCtrl()

Constant No  
External Access: Read/Write

*expOpenCloseCtrl - ISR\_ValveCtrl/Main - #100*

*expOpenCloseCtrl[0] - initVent/Main - \*#7*

*expOpenCloseCtrl[0] - ISR\_ValveCtrl/Main - \*#223*

expOpenCloseCtrlNew	INT[1000]	test
---------------------	-----------	------

Constant No  
External Access: Read/Write

*expOpenCloseCtrlNew - calcBreath/Main - \*#56, \*#74*

*expOpenCloseCtrlNew[0] - calcBreath/Main - #82*

expStartI	0	INT	test
-----------	---	-----	------

Constant No  
External Access: Read/Write

## LogixProgram\_L27\_Ver4.ACD

**expStartI (Continued)**

*expStartI - measArray/FirstScan - \*#36*  
*expStartI - measCalc/Main - #107, #23, #37, #53, #58, #78, #86, #97, #98, \*#30, \*#32, \*#54, \*#7*

**expTimeInMS**                    3000.0002                    REAL                    test

Constant                    No  
External Access:            Read/Write  
*expTimeInMS - calcBreath/Main - \*#15*

**F1\_pushed\_new**                0                    BOOL                    test

Constant                    No  
External Access:            Read/Write  
*F1\_pushed\_new - GUI/Handle\_GUI - #40, #42, #82*

**F1\_pushed\_old**                0                    BOOL                    test

Constant                    No  
External Access:            Read/Write  
*F1\_pushed\_old - GUI/Handle\_GUI - #40, #82, \*#41, \*#83*

**F2\_pushed\_new**                0                    BOOL                    test

Constant                    No  
External Access:            Read/Write  
*F2\_pushed\_new - GUI/Handle\_GUI - #129, #89, #91*

**F2\_pushed\_old**                0                    BOOL                    test

Constant                    No  
External Access:            Read/Write  
*F2\_pushed\_old - GUI/Handle\_GUI - #129, #89, \*#130, \*#90*

**F3\_pushed\_new**                0                    BOOL                    test

Constant                    No  
External Access:            Read/Write  
*F3\_pushed\_new - GUI/Handle\_GUI - #136, #138, #176*

**F3\_pushed\_old**                0                    BOOL                    test

Constant                    No  
External Access:            Read/Write  
*F3\_pushed\_old - GUI/Handle\_GUI - #136, #176, \*#137, \*#177*

**Fill**                            AOI\_FillArrayReal                    test

Constant                    No  
External Access:            Read/Write

*Fill - calcBreath/Main - \*#54, \*#55*

**Fill.EnableIn**                1                    BOOL

Enable Input - System Defined Parameter

**Fill.EnableOut**                1                    BOOL

Enable Output - System Defined Parameter

**Fill50**                            AOI\_FillArrayReal50                    test

Constant                    No  
External Access:            Read/Write

*Fill50 - initVent/initISR\_ValveCtrl - \*#5*

*Fill50 - measArray/FirstScan - \*#42, \*#43, \*#44, \*#45, \*#46, \*#47, \*#48, \*#49, \*#50, \*#51, \*#52*

**Fill50.EnableIn**                1                    BOOL

Enable Input - System Defined Parameter

**Fill50.EnableOut**                1                    BOOL

Enable Output - System Defined Parameter

**Fillint**                            AOI\_FillArrayInt                    test

Constant                    No  
External Access:            Read/Write

*Fillint - calcBreath/Main - \*#56, \*#74*

**Fillint.EnableIn**                1                    BOOL

Enable Input - System Defined Parameter

LogixProgram\_L27\_Ver4.ACD

**Fillint (Continued)**

**Fillint.EnableOut** 1 BOOL  
Enable Output - System Defined Parameter

**FillLoop** AOI\_FillArray\_Real2D test

Constant No  
External Access: Read/Write  
*FillLoop - MainFlow/Main - \*#56, \*#57, \*#58, \*#59, \*#60, \*#61, \*#62, \*#63, \*#64*  
*FillLoop - measArray/FirstScan - \*#8, \*#9*  
*FillLoop - measCalc/Main - \*#228, \*#229, \*#230*

**FillLoop.EnableIn** 1 BOOL  
Enable Input - System Defined Parameter

**FillLoop.EnableOut** 1 BOOL  
Enable Output - System Defined Parameter

**FiO2** INT test

25-100%  
Constant No  
External Access: Read/Write  
*FiO2 - calcBreath/Main - #36*  
*FiO2 - GUI/Handle\_GUI - #145, #146, #98, #99, \*#100, \*#145, \*#147, \*#98*  
*FiO2 - MainFlow/Main - #14, #9*  
*FiO2 - measArray/FirstScan - #61, \*#62, \*#64*

**FiO2Running** INT test

Constant No  
External Access: Read/Write  
*FiO2Running - MainFlow/Main - \*#14*  
*FiO2Running - measArray/FirstScan - #64*

**Floor** AOI\_Floor test

Constant No  
External Access: Read/Write  
*Floor - calcBreath/Main - \*#14, \*#17, \*#44, \*#49, \*#51*  
**Floor.EnableIn** 1 BOOL  
Enable Input - System Defined Parameter  
**Floor.EnableOut** 1 BOOL  
Enable Output - System Defined Parameter

**Flow** REAL[1000,2] measArray

Constant No  
External Access: Read/Write  
*Flow - measArray/measArrays - \*#32, \*#33*

**FlowAir** REAL[1000,2] test

1000 length 2D FIFO array with flow values and status bits indicating if values have been analysed (1) or not (0) by measCalc()  
Constant No  
External Access: Read/Write  
*FlowAir - ISR\_ValveCtrl/Main - #135, #142, #150, #163, #170, #179, #190, #198*  
*FlowAir - MainFlow/Main - \*#56*  
*FlowAir - measArray/FirstScan - \*#8*  
*FlowAir - measArray/measArrays - \*#34, \*#35*

**FlowArrMem** REAL[1000,2] test

Constant No  
External Access: Read/Write  
*FlowArrMem - ISR\_ValveCtrl/Main - \*#142, \*#163, \*#170, \*#179, \*#198*  
*FlowArrMem - MainFlow/Main - \*#59*  
*FlowArrMem - measCalc/Main - #108, #12, #18, #18, #27, #27, #41, #41, \*#228, \*#86*  
*FlowArrMem[0,0] - ISR\_ValveCtrl/Main - \*#135, \*#150, \*#190*

**flowCalcOn** BOOL test

Constant No  
External Access: Read/Write

LogixProgram\_L27\_Ver4.ACD

**flowCalcOn (Continued)**

*flowCalcOn - initVent/initISR\_ValveCtrl - \*#6  
flowCalcOn - ISR\_ValveCtrl/Main - #117*

<b>FlowImax</b>	18987.0	REAL	test
max insp flow in l/min			
Constant	No		
External Access:	Read/Write		
<i>FlowImax - measCalc/Main - #121, #165, *#87</i>			
<b>FlowImaxArr</b>		REAL[50]	test
50 length array with max insp flow in l/min, from measCalc()			
Constant	No		
External Access:	Read/Write		
<i>FlowImaxArr - measArray/FirstScan - *#44</i>			
<i>FlowImaxArr - measCalc/Main - #151, *#151, *#210</i>			
<i>FlowImaxArr[0] - measCalc/Main - *#165</i>			
<b>FlowImaxValue</b>		AOI_Max	test
Constant	No		
External Access:	Read/Write		
<i>FlowImaxValue - measCalc/Main - *#86</i>			
<b>FlowImaxValue.EnableIn</b>	0	BOOL	
Enable Input - System Defined Parameter			
<b>FlowImaxValue.EnableOut</b>	0	BOOL	
Enable Output - System Defined Parameter			
<b>FlowImaxValue.Out_Max</b>	18987.0	REAL	
<i>FlowImaxValue.Out_Max - measCalc/Main - #87</i>			
<b>FlowO2</b>		REAL[1000,2]	test
Constant	No		
External Access:	Read/Write		
<i>FlowO2 - ISR_ValveCtrl/Main - #136, #143, #151, #164, #171, #180, #191, #199</i>			
<i>FlowO2 - MainFlow/Main - *#57</i>			
<i>FlowO2 - measArray/measArrays - *#36, *#37</i>			
<b>FlowO2Mem</b>		REAL[1000,2]	test
Constant	No		
External Access:	Read/Write		
<i>FlowO2Mem - ISR_ValveCtrl/Main - *#143, *#164, *#171, *#180, *#199</i>			
<i>FlowO2Mem - MainFlow/Main - *#60</i>			
<i>FlowO2Mem - measCalc/Main - *#229</i>			
<i>FlowO2Mem[0,0] - ISR_ValveCtrl/Main - *#136, *#151, *#191</i>			
<b>FractionAir</b>	0.94936705	REAL	test
Constant	No		
External Access:	Read/Write		
<i>FractionAir - calcBreath/Main - #37, #41, *#36</i>			
<b>FractionO2</b>	0.050632954	REAL	test
Constant	No		
External Access:	Read/Write		
<i>FractionO2 - calcBreath/Main - #40, *#37</i>			
<b>GetSlaveValue</b>		MOTION_INSTRUCTION	test
Constant	No		
External Access:	Read/Write		
<i>GetSlaveValue - ISR_ValveCtrl/Main - *#61, *#63</i>			
<b>HMI_ScreenNo</b>	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>HMI_ScreenNo - GUI/Handle_GUI - *#12, *#3, *#7</i>			

LogixProgram\_L27\_Ver4.ACD

<b>HMIStart</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>HMIStart - GUI/Handle_GUI - #5</i>			
<b>HMIStop</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>HMIStop - GUI/Handle_GUI - #10</i>			
<b>i</b>	101	DINT	test
Constant	No		
External Access:	Read/Write		
<i>i - calcBreath/Main - #60, #61, #61, #62, #62, #63, #63, #64, #64, #65, #66, #67, #67, #68, #68, #69, *#59</i>			
<i>i - measCalc/Main - #189</i>			
<b>I_E</b>	'1:3'	STRING	test
part of I:E ratio: 1, 2, 3 or 4			
Constant	No		
External Access:	Read/Write		
<i>I_E - GUI/Handle_GUI - *#123, *#170</i>			
<b>I_E.LEN</b>	3	DINT	
part of I:E ratio: 1, 2, 3 or 4			
<b>I_E.DATA</b>		SINT	
part of I:E ratio: 1, 2, 3 or 4			
<b>IE</b>		STRING[4]	test
Constant	No		
External Access:	Read/Write		
<i>IE - GUI/Handle_GUI - #123, #170</i>			
<b>IE_counter</b>	2	DINT	test
Constant	No		
External Access:	Read/Write		
<i>IE_counter - GUI/Handle_GUI - #116, #117, #123, #163, #164, #170, *#116, *#118, *#163, *#165</i>			
<i>IE_counter - MainFlow/Main - #6</i>			
<b>IE_E</b>	3	INT	test
Constant	No		
External Access:	Read/Write		
<i>IE_E - calcBreath/Main - #11</i>			
<i>IE_E - MainFlow/Main - #18, #9, *#6</i>			
<i>IE_E - measArray/FirstScan - #85, *#86, *#88</i>			
<b>IE_ERunning</b>	3	INT	test
Constant	No		
External Access:	Read/Write		
<i>IE_ERunning - MainFlow/Main - *#18</i>			
<i>IE_ERunning - measArray/FirstScan - #88</i>			
<b>IE_I</b>	1	INT	test
Constant	No		
External Access:	Read/Write		
<i>IE_I - calcBreath/Main - #11, #11</i>			
<i>IE_I - MainFlow/Main - #17, #9</i>			
<i>IE_I - measArray/FirstScan - #79, *#80, *#82</i>			
<b>IE_IRunning</b>	1	INT	test
Constant	No		
External Access:	Read/Write		
<i>IE_IRunning - MainFlow/Main - *#17</i>			
<i>IE_IRunning - measArray/FirstScan - #82</i>			
<b>Insp</b>	1	BOOL	test

LogixProgram\_L27\_Ver4.ACD

Insp (Continued)

Constant No  
External Access: Read/Write  
*Insp - ISR\_ValveCtrl/Main - #26, #6, #6, \*#114, \*#12, \*#16, \*#71, \*#9*  
*Insp - MainFlow/Main - \*#66*

inspAirFlowCtrl

REAL[650]

test

650 length array with flow values for insp air valve, moved from nextInspAirFlowCtrl[] by ISR\_ValveCtrl()  
Constant No  
External Access: Read/Write  
*inspAirFlowCtrl - ISR\_ValveCtrl/Main - #63*  
*inspAirFlowCtrl[0] - initVent/Main - \*#6*  
*inspAirFlowCtrl[0] - ISR\_ValveCtrl/Main - \*#87*

inspAirFlowCtrlNew

REAL[650]

test

Constant No  
External Access: Read/Write  
*inspAirFlowCtrlNew - calcBreath/Main - \*#55, \*#63, \*#66, \*#69*  
*inspAirFlowCtrlNew[0] - calcBreath/Main - #78*

inspCount

0

INT

test

between 1-650 to indicate where valveCtrl should read from inspAirFlowCtrl[] and inspO2FlowCtrl[]  
Constant No  
External Access: Read/Write  
*inspCount - ISR\_ValveCtrl/Main - #11, #41, #61, #63, #7, #70, \*#41, \*#97*

inspCountLength

100

INT

test

between 1-650 to indicate when valveCtrl should stop reading from inspAirFlowCtrl[] and inspO2FlowCtrl[], copied from nextInspCountLength by ISR\_ValveCtrl()

Constant No  
External Access: Read/Write

*inspCountLength - initVent/Main - #11, \*#8*  
*inspCountLength - ISR\_ValveCtrl/Main - #11, #7, #70, \*#88*

inspCountLengthNew

100.0

REAL

test

Constant No  
External Access: Read/Write

*inspCountLengthNew - calcBreath/Main - #16, #59, #79, \*#14*

inspCounts

100.00001

REAL

test

Constant No  
External Access: Read/Write

*inspCounts - calcBreath/Main - #14, \*#13*

InspiratoryRampTimeInMS

200.00002

REAL

test

Constant No  
External Access: Read/Write

*InspiratoryRampTimeInMS - calcBreath/Main - #22, #27, #43, \*#21*

inspO2FlowCtrl

REAL[650]

test

650 length array with flow values for insp O2 valve, moved from nextInspO2FlowCtrl[] by ISR\_ValveCtrl()  
Constant No  
External Access: Read/Write  
*inspO2FlowCtrl - ISR\_ValveCtrl/Main - #61*  
*inspO2FlowCtrl[0] - initVent/Main - \*#5*  
*inspO2FlowCtrl[0] - ISR\_ValveCtrl/Main - \*#86*

inspO2FlowCtrlNew

REAL[650]

test

Constant No  
External Access: Read/Write

*inspO2FlowCtrlNew - calcBreath/Main - \*#54, \*#62, \*#65, \*#68*  
*inspO2FlowCtrlNew[0] - calcBreath/Main - #77*

## LogixProgram\_L27\_Ver4.ACD

■ <b>inspPauseCalc</b>	40.000004	REAL	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseCalc - calcBreath/Main</i> - #51, *#50			
■ <b>inspPauseNrofAdjustments</b>	40.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseNrofAdjustments - calcBreath/Main</i> - #67, *#51			
■ <b>inspPauseStartI</b>	0	INT	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseStartI - measArray/FirstScan</i> - *#35			
<i>inspPauseStartI - measCalc/Main</i> - #120, #130, #94, *#29			
■ <b>inspPauseTime</b>	10	INT	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseTime - calcBreath/Main</i> - #8			
<i>inspPauseTime - MainFlow/Main</i> - #19, #9			
<i>inspPauseTime - measArray/FirstScan</i> - #91, *#92, *#94			
■ <b>inspPauseTimeInMS</b>	400.00003	REAL	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseTimeInMS - calcBreath/Main</i> - #22, #50, *#20			
<i>inspPauseTimeInMS - measCalc/Main</i> - #95			
■ <b>inspPauseTimePercent</b>	0.1	REAL	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseTimePercent - calcBreath/Main</i> - #20, *#8			
■ <b>inspPauseTimeRunning</b>	10	INT	test
Constant	No		
External Access:	Read/Write		
<i>inspPauseTimeRunning - MainFlow/Main</i> - *#19			
<i>inspPauseTimeRunning - measArray/FirstScan</i> - #94			
■ <b>inspRampTime</b>	5	INT	test
Constant	No		
External Access:	Read/Write		
<i>inspRampTime - calcBreath/Main</i> - #9			
<i>inspRampTime - MainFlow/Main</i> - #20, #9			
<i>inspRampTime - measArray/FirstScan</i> - #97, *#100, *#98			
■ <b>inspRampTimePercent</b>	0.05	REAL	test
Constant	No		
External Access:	Read/Write		
<i>inspRampTimePercent - calcBreath/Main</i> - #21, *#9			
■ <b>inspRampTimeRunning</b>	5	INT	test
Constant	No		
External Access:	Read/Write		
<i>inspRampTimeRunning - MainFlow/Main</i> - *#20			
<i>inspRampTimeRunning - measArray/FirstScan</i> - #100			
■ <b>inspStartI</b>	0	INT	test
Constant	No		
External Access:	Read/Write		
<i>inspStartI - measArray/FirstScan</i> - *#34			
<i>inspStartI - measCalc/Main</i> - #107, #139, #14, #23, #37, #49, #54, #64, #68, #78, #86, *#19, *#50, *#6			

LogixProgram\_L27\_Ver4.ACD

■ <b>inspTimeinMS</b>	1000.00006	REAL	test
Constant	No		
External Access:	Read/Write		
<i>inspTimeinMS - calcBreath/Main</i> - #13, #15, #22, *#11			
■ <b>last_state</b>	3	DINT	test
Constant	No		
External Access:	Read/Write		
<i>last_state - GUI/Handle_GUI</i> - #50, *#186			
■ <b>Local:1:O</b>		AB:Embedded_DiscreteIO1:O:0	test
Constant	No		
External Access:	Read/Write		
<b>Local:1:O.Data.0</b>	1	BOOL	
<i>expOpen - ISR_ValveCtrl/Main</i> - *#100, *#54, *#60			
■ <b>Local:2:I</b>		AB:Embedded_AnalogIO1:I:0	test
Constant	No		
External Access:	Read/Write		
<b>Local:2:I.Ch0Data</b>	2#0000_0011_1101_0101	INT	
<i>AnalogInputFlowO2 - measArray/measArrays</i> - #111, #36			
<b>Local:2:I.Ch1Data</b>	2#0000_0011_1101_0100	INT	
<i>AnalogInputFlowAir - measArray/measArrays</i> - #109, #34			
<b>Local:2:I.Ch2Data</b>	2#0000_1111_0110_0111	INT	
<i>AnalogInputPressure - ISR_ValveCtrl/Main</i> - #1			
<i>AnalogInputPressure - measArray/measArrays</i> - #110, #38			
<b>Local:2:I.Ch3Data</b>	2#0000_0100_0010_0011	INT	
<i>AnalogInputFlow - measArray/measArrays</i> - #32			
■ <b>Local:2:O</b>		AB:Embedded_AnalogIO1:O:0	test
Constant	No		
External Access:	Read/Write		
<b>Local:2:O.Ch0Data</b>	0	INT	
<i>AnalogOutFlowO2 - ISR_ValveCtrl/Main</i> - *#52, *#62, *#98			
<b>Local:2:O.Ch1Data</b>	0	INT	
<i>AnalogOutFlowAir - ISR_ValveCtrl/Main</i> - *#53, *#64, *#99			
■ <b>maxFlowAllowed</b>	100	INT	test
Constant	No		
External Access:	Read/Write		
<i>maxFlowAllowed - calcBreath/Main</i> - #31			
<i>maxFlowAllowed - measArray/FirstScan</i> - *#58			
■ <b>Mean50</b>		AOI_Mean50	test
Constant	No		
External Access:	Read/Write		
<i>Mean50 - measCalc/Main</i> - *#206, *#208, *#210, *#212, *#214, *#216, *#218, *#220, *#222, *#224			
<b>Mean50.EnableIn</b>	0	BOOL	
Enable Input - System Defined Parameter			
<b>Mean50.EnableOut</b>	0	BOOL	
Enable Output - System Defined Parameter			
<b>Mean50.Out_Mean</b>	2.059567	REAL	
<i>Mean50.Out_Mean - measCalc/Main</i> - #207, #209, #211, #213, #215, #217, #219, #221, #223, #225			
■ <b>message</b>	'Data Ready'	STRING	test
Constant	No		
External Access:	Read/Write		
<i>message - calcBreath/Main</i> - *#32			
<i>message - GUI/Handle_GUI</i> - *#184, *#196			
<i>message - initVent/Main</i> - *#16			
■ <b>Messages_to_user</b>		STRING[10]	test
Constant	No		
External Access:	Read/Write		

LogixProgram\_L27\_Ver4.ACD**Messages\_to\_user (Continued)**

Messages\_to\_user[0] - GUI/Handle\_GUI - #184  
Messages\_to\_user[1] - initVent/Main - #16  
Messages\_to\_user[2] - calcBreath/Main - #32

<b>Mode</b>	'VC'	STRING	test
Constant	No		
External Access:	Read/Write		
Mode - measArray/FirstScan - *#33			
Mode - measCalc/Main - #28			
<b>MV</b>	48025.97	REAL	test
minute ventilation in l/min			
Constant	No		
External Access:	Read/Write		
MV - measCalc/Main - #169, *#113			
<b>MVArr</b>		REAL[50]	test
50 length array with minute ventilation in l/min, from measCalc()			
Constant	No		
External Access:	Read/Write		
MVArr - measArray/FirstScan - *#48			
MVArr - measCalc/Main - #155, *#155, *#218			
MVArr[0] - measCalc/Main - *#169			
<b>NewExpOpenCloseAvail</b>	0	BOOL	test
signal to valveCtrl to include new expOpenCloseCtrl[] to overwrite when expiration ended, set by calcBreath()			
Constant	No		
External Access:	Read/Write		
NewExpOpenCloseAvail - calcBreath/Main - *#87			
NewExpOpenCloseAvail - ISR_ValveCtrl/Main - #215, *#225			
<b>newInspFlowAvil</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
newInspFlowAvil - calcBreath/Main - *#86			
newInspFlowAvil - ISR_ValveCtrl/Main - #74, *#89			
<b>NewSettings</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
NewSettings - calcBreath/Main - *#85			
NewSettings - MainFlow/Main - #42, *#13			
NewSettings - measArray/FirstScan - *#16			
<b>NewTimeRate_p5</b>	10000	DINT	test
Constant	No		
External Access:	Read/Write		
NewTimeRate_p5 - measArray/measArrays - #22, *#19			
<b>NewTimeRate_p8</b>	10000	DINT	test
Constant	No		
External Access:	Read/Write		
NewTimeRate_p8 - measArray/measArrays - #22, *#20			
<b>nextExpCountLength</b>	300	INT	test
between 1-1000 to indicate when valveCtrl should stop reading from expOpenCloseCtrl[], set by calcBreath()			
Constant	No		
External Access:	Read/Write		
nextExpCountLength - calcBreath/Main - *#83			
nextExpCountLength - initVent/Main - #9			
nextExpCountLength - ISR_ValveCtrl/Main - #224			
<b>nextExpOpenCloseCtrl</b>		INT[1000]	test

## LogixProgram\_L27\_Ver4.ACD

**nextExpOpenCloseCtrl (Continued)**

1000 length array with open (1) or close (0) for exp valve, set by calcBreath()  
Constant No  
External Access: Read/Write  
*nextExpOpenCloseCtrl[0]* - calcBreath/Main - \*#82  
*nextExpOpenCloseCtrl[0]* - initVent/Main - #7  
*nextExpOpenCloseCtrl[0]* - ISR\_ValveCtrl/Main - #223

**nextInspAirFlowCtrl**

REAL[650] test

650 length array with flow values for insp air valve, set by calcBreath()  
Constant No  
External Access: Read/Write  
*nextInspAirFlowCtrl[0]* - calcBreath/Main - \*#78  
*nextInspAirFlowCtrl[0]* - initVent/Main - #6  
*nextInspAirFlowCtrl[0]* - ISR\_ValveCtrl/Main - #87

**nextInspCountLength**100 INT test  
between 1-650 to indicate when valveCtrl should stop reading from inspAirFlowCtrl[] and inspO2FlowCtrl[], set by calcBreath()

Constant No  
External Access: Read/Write  
*nextInspCountLength* - calcBreath/Main - \*#79  
*nextInspCountLength* - initVent/Main - #8  
*nextInspCountLength* - ISR\_ValveCtrl/Main - #88

**nextInspO2FlowCtrl**

REAL[650] test

650 length array with flow values for insp O2 valve, set by calcBreath()  
Constant No  
External Access: Read/Write  
*nextInspO2FlowCtrl[0]* - calcBreath/Main - \*#77  
*nextInspO2FlowCtrl[0]* - initVent/Main - #5  
*nextInspO2FlowCtrl[0]* - ISR\_ValveCtrl/Main - #86

**O2ConstantFlow**

3.6455727 REAL test

Constant No  
External Access: Read/Write  
*O2ConstantFlow* - calcBreath/Main - #45, #65, \*#40

**O2Down**

CAM\_PROFILE[1010] test

Constant No  
External Access: Read/Write  
*O2Down* - ValveProfile/Main - 0(MAPC)

**O2mAder**

0.0 REAL test

Constant No  
External Access: Read/Write  
*O2mAder* - ISR\_ValveCtrl/Main - \*#61

**O2mAOut**

0.0 REAL test

Constant No  
External Access: Read/Write  
*O2mAOut* - ISR\_ValveCtrl/Main - #62, \*#61

**O2mASlope**

416748.25 REAL test

Constant No  
External Access: Read/Write  
*O2mASlope* - ISR\_ValveCtrl/Main - \*#61

**O2RampFlowIncreasePerAdjust**

0.18227863 REAL test

Constant No  
External Access: Read/Write  
*O2RampFlowIncreasePerAdjust* - calcBreath/Main - #62, \*#45

**O2Up**

CAM\_PROFILE[1010] test

Constant No

LogixProgram\_L27\_Ver4.ACD

**O2Up (Continued)**

External Access: Read/Write  
*O2Up - ISR\_ValveCtrl/Main - #61*  
*O2Up - ValveProfile/Main - 0(MAPC)*

**PEEP** 4402.0 REAL test  
positive end expiratory pressure in cmH2O  
Constant No  
External Access: Read/Write  
*PEEP - measCalc/Main - #127, #131, #164, \*#68*

**PEEPArr** REAL[50] test  
50 length array with positive end expiratory pressure in cmH2O, from measCalc()  
Constant No  
External Access: Read/Write  
*PEEPArr - measArray/FirstScan - \*#43*  
*PEEPArr - measCalc/Main - #150, \*#150, \*#208*  
*PEEPArr[0] - measCalc/Main - \*#164*

**Ppeak** 10534.0 REAL test  
Ppeak  
Constant No  
External Access: Read/Write  
*Ppeak - measCalc/Main - #121, #127, #163, \*#79*

**PpeakArr** REAL[50] test  
50 length array with Ppeak from measCalc(), used for minAvg  
Constant No  
External Access: Read/Write  
*PpeakArr - measArray/FirstScan - \*#42*  
*PpeakArr - measCalc/Main - #149, \*#149, \*#206*  
*PpeakArr[0] - measCalc/Main - \*#163*

**Ppeak.MaxValue** AOI\_Max test  
Constant No  
External Access: Read/Write  
*Ppeak.MaxValue - measCalc/Main - \*#78*  
**Ppeak.MaxValue.EnableIn** 0 BOOL  
Enable Input - System Defined Parameter  
**Ppeak.MaxValue.EnableOut** 0 BOOL  
Enable Output - System Defined Parameter  
**Ppeak.MaxValue.Out\_Max** 10534.0 REAL  
*Ppeak.MaxValue.Out\_Max - measCalc/Main - #79*

**Pplat** 9761.556 REAL test  
plateau insp pressure in cmH2O  
Constant No  
External Access: Read/Write  
*Pplat - measCalc/Main - #121, #131, #166, \*#100, \*#102*

**PplatArr** REAL[50] test  
50 length array with plateau insp pressure in cmH2O, from measCalc()  
Constant No  
External Access: Read/Write  
*PplatArr - measArray/FirstScan - \*#45*  
*PplatArr - measCalc/Main - #152, \*#152, \*#212*  
*PplatArr[0] - measCalc/Main - \*#166*

**Pplat.BreathFrac** 0.2 REAL test  
Constant No  
External Access: Read/Write  
*Pplat.BreathFrac - measArray/FirstScan - \*#40*  
*Pplat.BreathFrac - measCalc/Main - #95*

LogixProgram\_L27\_Ver4.ACD

<b>PplatMean</b>		AOI_Mean2D	test
Constant	No		
External Access:	Read/Write		
<i>PplatMean - measCalc/Main - *#99</i>			
<b>PplatMean.EnableIn</b>	0	BOOL	
Enable Input - System Defined Parameter			
<b>PplatMean.EnableOut</b>	0	BOOL	
Enable Output - System Defined Parameter			
<b>PplatMean.Out_Mean</b>	9761.556	REAL	
<i>PplatMean.Out_Mean - measCalc/Main - #100</i>			
<b>PressArr</b>		REAL[1000,2]	test
1000 length 2D FIFO array with pressure values and status bits indicating if values have been analysed (1) or not (0) by measCalc()			
Constant	No		
External Access:	Read/Write		
<i>PressArr - ISR_ValveCtrl/Main - #137, #144, #152, #165, #172, #181, #192, #200</i>			
<i>PressArr - MainFlow/Main - *#58</i>			
<i>PressArr - measArray/FirstScan - *#9</i>			
<i>PressArr - measArray/measArrays - *#38, *#39</i>			
<b>PressArrMem</b>		REAL[1000,2]	test
Constant	No		
External Access:	Read/Write		
<i>PressArrMem - ISR_ValveCtrl/Main - *#144, *#165, *#172, *#181, *#200</i>			
<i>PressArrMem - MainFlow/Main - *#61</i>			
<i>PressArrMem - measCalc/Main - #68, *#230, *#78, *#99</i>			
<i>PressArrMem[0,0] - ISR_ValveCtrl/Main - *#137, *#152, *#192</i>			
<b>PressHthresholdUser1</b>	35.0	REAL	ISR_ValveCtrl
Constant	No		
External Access:	Read/Write		
<i>PressHthresholdUser1 - ISR_ValveCtrl/Main - #43</i>			
<b>Pressure</b>	-0.30625	REAL	test
Constant	No		
External Access:	Read/Write		
<i>Pressure - ISR_ValveCtrl/Main - #43, *#1</i>			
<b>RampNrCalc</b>	20.000002	REAL	test
Constant	No		
External Access:	Read/Write		
<i>RampNrCalc - calcBreath/Main - #44, *#43</i>			
<b>RampNrOfAdjustments</b>	20.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>RampNrOfAdjustments - calcBreath/Main - #45, #46, #60, #64, #64, #67, #67, *#44</i>			
<b>RampTime</b>		REAL[650]	test
Constant	No		
External Access:	Read/Write		
<i>RampTime - calcBreath/Main - *#61</i>			
<b>Raw</b>	2.440968	REAL	test
airway resistance in cmH2O/l/s			
Constant	No		
External Access:	Read/Write		
<i>Raw - measCalc/Main - #172, *#121, *#123</i>			
<b>RawArr</b>		REAL[50]	test
50 length array with airway resistance in cmH2O/l/s			
Constant	No		
External Access:	Read/Write		
<i>RawArr - measArray/FirstScan - *#51</i>			

## LogixProgram\_L27\_Ver4.ACD

**RawArr (Continued)**

RawArr - measCalc/Main - #158, \*#158, \*#224  
RawArr[0] - measCalc/Main - \*#172

**Round** AOI\_Round

test

Constant No  
External Access: Read/Write

**Round.EnableIn** 0 BOOL

Enable Input - System Defined Parameter

**Round.EnableOut** 0 BOOL

Enable Output - System Defined Parameter

**RR** REAL

test

respiratory rate in 1/min  
Constant No  
External Access: Read/Write

RR - measCalc/Main - #113, #167, \*#65

**RRArr** REAL[50]

test

50 length array with respiratory rate in 1/min, from measCalc()

Constant No  
External Access: Read/Write

RRArr - measArray/FirstScan - \*#46

RRArr - measCalc/Main - #153, \*#153, \*#214

RRArr[0] - measCalc/Main - \*#167

**RRate** INT

test

Constant No  
External Access: Read/Write

RRate - calcBreath/Main - #10

RRate - GUI/Handle\_GUI - #110, #111, #157, #158, \*#110, \*#112, \*#157, \*#159

RRate - MainFlow/Main - #16, #9

RRate - measArray/FirstScan - #73, \*#74, \*#76

**RRateRunning** INT

test

Constant No  
External Access: Read/Write

RRateRunning - MainFlow/Main - \*#16

RRateRunning - measArray/FirstScan - #76

**RunModes** STRING[4]

test

Constant No  
External Access: Read/Write

RunModes[0] - measArray/FirstScan - #33

RunModes[0] - measCalc/Main - #28

**SampelRateChangeAllowed** 0 BOOL

test

Constant No  
External Access: Read/Write

SampelRateChangeAllowed - measArray/FirstScan - \*#4

SampelRateChangeAllowed - measArray/measArrays - \*#23

**Samplerate** DINT

test

The rate in Hz that samples is collected

Constant No  
External Access: Read/Write

Samplerate - measArray/FirstScan - \*#5

Samplerate - measArray/measArrays - #12, #15, #18, #6, #8, \*#13, \*#9

Samplerate - measCalc/Main - #108, #64, #95

**Samplerate\_Old** DINT

test

Constant No  
External Access: Read/Write

## LogixProgram\_L27\_Ver4.ACD

**Samplerate\_Old (Continued)**

*Samplerate\_Old - measArray/FirstScan - \*#6*  
*Samplerate\_Old - measArray/measArrays - #6, \*#18, \*#24*

SetupDone	1	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>SetupDone - MainFlow/Main - #29</i>			
<i>SetupDone - measArray/FirstScan - *#15</i>			
<i>SetupDone - Setup/Main - *#10</i>			
SetupState	0	DINT	test
Constant	No		
External Access:	Read/Write		
<i>SetupState - MainFlow/Main - #35, *#31, *#36</i>			
<i>SetupState - measArray/FirstScan - *#23</i>			
<i>SetupState - Setup/Main - #1, *#11, *#4, *#7</i>			
Start_Stop	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>Start_Stop - GUI/Handle_GUI - #16, *#11, *#6</i>			
StartPoint	186	INT	test
Constant	No		
External Access:	Read/Write		
<i>StartPoint - measCalc/Main - #99, *#97</i>			
StartPointMean	8.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>StartPointMean - measCalc/Main - #97, *#96</i>			
StartPointValue	8.2	REAL	test
Constant	No		
External Access:	Read/Write		
<i>StartPointValue - measCalc/Main - #96, *#95</i>			
StartSendPosition	-30	DINT	test
Constant	No		
External Access:	Read/Write		
<i>StartSendPosition - ISR_ValveCtrl/Main - #131, #133, #134, #141, #148, #188, #189, #196, *#129, *#187</i>			
StartVent	0	BOOL	test
1: start ventilating patient, 0: do nothing			
Constant	No		
External Access:	Read/Write		
<i>StartVent - GUI/Handle_GUI - *#19</i>			
<i>StartVent - initVent/Main - *#14, *#25</i>			
<i>StartVent - MainFlow/Main - #78</i>			
<i>StartVent - measArray/FirstScan - *#17</i>			
state	3	DINT	test
Constant	No		
External Access:	Read/Write		
<i>state - GUI/Handle_GUI - #144, #186, #53, #54, #59, #97, *#50, *#54, *#55</i>			
state_selected	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>state_selected - GUI/Handle_GUI - #138, #49, #91, *#187, *#51</i>			
state0	0	BOOL	test
Constant	No		

## LogixProgram\_L27\_Ver4.ACD

**state0 (Continued)**

External Access: Read/Write  
*state0 - GUI/Handle\_GUI - \*#188, \*#60, \*#64, \*#68, \*#72*

**state1** 0 BOOL test

Constant No  
External Access: Read/Write  
*state1 - GUI/Handle\_GUI - \*#189, \*#61, \*#65, \*#69, \*#73*

**state2** 0 BOOL test

Constant No  
External Access: Read/Write  
*state2 - GUI/Handle\_GUI - \*#190, \*#62, \*#66, \*#70, \*#74*

**state3** 0 BOOL test

Constant No  
External Access: Read/Write  
*state3 - GUI/Handle\_GUI - \*#191, \*#63, \*#67, \*#71, \*#75*

**StateCalcBreath** 0 DINT test

Constant No  
External Access: Read/Write  
*StateCalcBreath - calcBreath/Main - #5, \*#88*  
*StateCalcBreath - MainFlow/Main - #44, \*#45*  
*StateCalcBreath - measArray/FirstScan - \*#56*

**StateInitVent** 0 DINT test

Constant No  
External Access: Read/Write  
*StateInitVent - initVent/Main - #1, \*#12, \*#26*  
*StateInitVent - MainFlow/Main - #80, \*#81, \*#87*

**StateMain** 0 DINT test

Constant No  
External Access: Read/Write  
*StateMain - MainFlow/Main - #2, #26, \*#3, \*#30, \*#48, \*#69, \*#74, \*#86*  
*StateMain - measArray/FirstScan - \*#14*

**StateValveCtrl** 0 DINT test

Constant No  
External Access: Read/Write  
*StateValveCtrl - ISR\_ValveCtrl/Main - #113, #222, #23, #39, #50, #58, #68, #84, #94, \*#103, \*#109, \*#216, \*#218, \*#226, \*#27, \*#35, \*#44,*  
*\*#46, \*#55, \*#65, \*#75, \*#77, \*#80, \*#90*

**StopVent** 0 BOOL test

1: stop ventilating patient after finished breath, shut down not necessary ISRs  
Constant No  
External Access: Read/Write  
*StopVent - GUI/Handle\_GUI - \*#25*  
*StopVent - MainFlow/Main - #52, \*#65*  
*StopVent - measArray/FirstScan - \*#18*

**Time** 4800 DINT test

Constant No  
External Access: Read/Write  
*Time - GUI/Handle\_GUI - #196, \*#195*

**Timer** 0 DINT test

Constant No  
External Access: Read/Write  
*Timer - GUI/Handle\_GUI - #195, #32, #33, \*#141, \*#32, \*#36, \*#45, \*#94*

**TimeRate** 10000 DINT test

Constant No

## LogixProgram\_L27\_Ver4.ACD

**TimeRate (Continued)**

External Access: Read/Write  
*TimeRate - measArray/measArrays - #16, #17, #22, #22, \*#15*  
*TimeRate - measCalc/Main - #30*

■ <b>TimerDone</b>	0	BOOL	test
Constant	No		
External Access:	Read/Write		
<i>TimerDone - GUI/Handle_GUI - #183, *#192, *#197, *#35</i>			
■ <b>tipause</b>	250	INT	test
Constant	No		
External Access:	Read/Write		
<i>tipause - measCalc/Main - #30</i>			
■ <b>totalCount</b>	400.00003	REAL	test
Constant	No		
External Access:	Read/Write		
<i>totalCount - calcBreath/Main - #16, *#12</i>			
■ <b>totTime</b>	0	INT	test
Constant	No		
External Access:	Read/Write		
<i>totTime - measArray/FirstScan - *#38</i>			
<i>totTime - measCalc/Main - #195, #197, *#195</i>			
■ <b>valveUpdInterval</b>	10.0	REAL	test
Constant	No		
External Access:	Read/Write		
<i>valveUpdInterval - calcBreath/Main - #12, #13, #43, #48, #50, #61</i>			
<i>valveUpdInterval - measArray/FirstScan - *#57</i>			
■ <b>VentOn</b>	0	BOOL	test
1: keep ventilatin patient, checked by ISR_ValveCtrl() to stop after breath and measCalc()			
Constant	No		
External Access:	Read/Write		
<i>VentOn - GUI/Handle_GUI - #18, #24</i>			
<i>VentOn - initVent/Main - *#15, *#24</i>			
<i>VentOn - ISR_ValveCtrl/Main - #21</i>			
<i>VentOn - MainFlow/Main - *#55</i>			
<i>VentOn - measArray/FirstScan - *#19</i>			
■ <b>Vt</b>	600	INT	test
200-600 ml			
Constant	No		
External Access:	Read/Write		
<i>Vt - calcBreath/Main - #27</i>			
<i>Vt - GUI/Handle_GUI - #104, #105, #151, #152, *#104, *#106, *#151, *#153</i>			
<i>Vt - MainFlow/Main - #15, #9</i>			
<i>Vt - measArray/FirstScan - #67, *#68, *#70</i>			
■ <b>Vti</b>	3201.7312	REAL	test
insp tidal volume in ml			
Constant	No		
External Access:	Read/Write		
<i>Vti - measCalc/Main - #108, #113, #127, #131, #168, *#108</i>			
■ <b>VtiArr</b>		REAL[50]	test
50 length array with insp tidal volume in ml, from measCalc()			
Constant	No		
External Access:	Read/Write		
<i>VtiArr - measArray/FirstScan - *#47</i>			
<i>VtiArr - measCalc/Main - #154, *#154, *#216</i>			
<i>VtiArr[0] - measCalc/Main - *#168</i>			

LogixProgram\_L27\_Ver4.ACD

<b>VtRunning</b>	600	INT	test
Constant	No		
External Access:	Read/Write		
<i>VtRunning - MainFlow/Main - *#15</i>			
<i>VtRunning - measArray/FirstScan - #70</i>			

**General****Configuration**

Type:	Periodic	Watchdog:	300.000 ms
Period:	200.000 ms	Disable automatic output processing to reduce task overhead:	No
Priority:	10	Inhibit task:	No

**Program Schedule****Scheduled**

GUI ValveProfile

**Unscheduled****Monitor****Scan Times(Elapsed Time)**

Max: 0.486000 ms Last: 0.299000 ms

**Interval Times(Elapsed Times Between Triggers)**

Max: 200.207993 ms Min: 199.792007 ms

Task overlap count: 0

## General

### Configuration

Main:	 Handle_GUI	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

### Monitor

#### Scan Times(Execution Time)

Max:	67 us	Last:	36 us
------	-------	-------	-------

```
// Start and stop of ventilator

HMI_ScreenNo :=0;

If HMISstart then
  Start_Stop :=1;
  HMI_ScreenNo :=2;
end_if;

If HMISstop then
  Start_Stop :=0;
  HMI_ScreenNo :=1;
End_if;

if Start_Stop then                                // If momantary push button from GUI is on

  if not VentOn then;
    StartVent := 1;
  end_if;

  // If VentOn from ISR_ValveCtrl is NOT on
  // Set start signal to MainFlow sequence

else

  if VentOn then;
    StopVent :=1;
  end_if;
end_if;

// Control of the F1 push button - Shifts between the 4 input just above the Fx buttons
if data_changed then                                // start timer if data is changed
  Timer := Timer+1;
  if Timer >= 25 then;

    TimerDone := 1;
    Timer := 0;
  end_if;
end_if;

IF F1_pushed_new AND NOT F1_pushed_old Then
  F1_pushed_old := 1;
  IF F1_pushed_new THEN
    IF data_changed THEN
      data_changed := 0;
      Timer := 0;
    end_if;
    END_IF;
    data_changed :=1;
    IF NOT(state_selected) THEN
      state := last_state;
      state_selected := 1;
    ELSE
      IF state <3 THEN
        state := (state+1);
      ELSE state := 0;
    end_if;

    END_IF;
    CASE state OF
```

```
0: state0 := 1;
state1 := 0;
state2 := 0;
state3 := 0;
1: state0 := 0;
state1 := 1;
state2 := 0;
state3 := 0;
2: state0 := 0;
state1 := 0;
state2 := 1;
state3 := 0;
3: state0 := 0;
state1 := 0;
state2 := 0;
state3 := 1;
ELSE;
END_CASE;
END_IF;
END_IF;

if not F1_pushed_new AND F1_pushed_old then;
  F1_pushed_old :=0;
end_if;

// Control of the F2 button - Add value to selected variable

IF F2_pushed_new AND NOT F2_pushed_old Then
  F2_pushed_old := 1;
  IF (F2_pushed_new AND state_selected) THEN          // increment value if state is selected
    IF data_changed THEN                            // reset timer if data if changed repeatedly
      data_changed :=0;
      Timer := 0;                                  // reset timer
    END_IF;
    data_changed := 1;                            // signal the timer to start
    CASE state OF
      0: Fi02 := Fi02 + 5;                      // increment value in case the O2 is selected
      if Fi02 >= 100 then;
        Fi02 := 100;
      else;
        end_if;
      1: Vt := Vt + 50;                         // increment value in case the Vt is selected
      if Vt >= 600 then;
        Vt := 600;
      else;
        end_if;
      2: RRate := RRate + 1;                     // increment value in case the RRin is selected
      if RRate >= 35 then;
        RRate := 35;
      else;
        end_if;
      3: IE_counter := IE_counter + 1;           // increment value in case the I:E is selected
      if IE_counter >= 2 then;
        IE_counter := 2;
```

```
        else;

        end_if;
        COP(IE[IE_counter], I_E, 82);
    ELSE;
    END_CASE;
END_IF;
END_IF;

if not F2_pushed_new AND F2_pushed_old then;
    F2_pushed_old :=0;
end_if;

// Control of the F3 button - Sub value to selected variable

IF F3_pushed_new AND NOT F3_pushed_old Then
    F3_pushed_old := 1;
    IF (F3_pushed_new AND state_selected) THEN
        IF data_changed THEN
            data_changed := 0;
            Timer := 0;
        END_IF;
        data_changed := 1;
        CASE state OF
            0: Fi02 := Fi02 - 5;           // decrement value if state is selected
                if Fi02 <= 0 then;
                    Fi02 := 0;
                else;
                    //Fi02 := 0;
                end_if;
            1: Vt := Vt - 50;           // decrement value in case the Vt is selected
                if Vt <= 200 then;
                    Vt := 200;
                else;
                    end_if;
            2: RRate := RRate - 1;       // decrement value in case the RRin is selected
                if RRate <= 6 then;
                    RRate := 6;
                else;
                    end_if;
            3: IE_counter := IE_counter - 1; // increment value in case the I:E is selected
                if IE_counter <= 0 then;
                    IE_counter := 0;
                    //COP(IE[IE_counter], I_E, 82);
                else;
                    end_if;
                    COP(IE[IE_counter], I_E, 82);
    ELSE;
    END_CASE;
END_IF;
END_IF;

if      F3_pushed_new AND F3_pushed_old then;
nbt_pushed_old :=0;
```

```
end_if;

// If timer is done then deselect all states for new start

IF TimerDone THEN
    COP(Messages_to_user[0], message, 82);
    data_changed := 0;
    last_state := state;
    state_selected := 0;
    state0 := 0;
    state1 := 0;
    state2 := 0;
    state3 := 0;
    TimerDone := 0;
ELSE
    IF data_changed THEN
        Time := Timer*200;
        COP(Time, message, 82);
        TimerDone := 0;
    END_IF;
END_IF;
```

**General**

Type:  Structured Text      Number of Lines: Original 203  
In Program:  GUI

## General


### Configuration

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

## Monitor

### Scan Times(Execution Time)

Max:	77 us	Last:	50 us
------	-------	-------	-------





**General**

Type:  Ladder Diagram      Number of Rungs: 1  
In Program:  ValveProfile

**General**

Get raw analog input and control of the breathing - MUST HAVE HIGHEST priority of all peridoc tasks

**Configuration**

Type:	Periodic	Watchdog:	2000.000 ms
Period:	10.000 ms	Disable automatic output processing to reduce task overhead:	No
Priority:	5	Inhibit task:	No

**Program Schedule****Scheduled**

measArray	ISR_ValveCtrl
-----------	---------------

**Unscheduled****Monitor****Scan Times(Elapsed Time)**

Max:	1.412000 ms	Last:	1.258000 ms
------	-------------	-------	-------------

**Interval Times(Elapsed Times Between Triggers)**

Max:	10.943000 ms	Min:	9.867000 ms
------	--------------	------	-------------

Task overlap count:

0

## General

### Configuration

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

## Monitor

### Scan Times(Execution Time)

Max:	814 us	Last:	420 us
------	--------	-------	--------

```
// From Analog_values

counter :=0;
SampelRateChangeAllowed := 0;
Samplerate := 100 ;
Samplerate_Old := 0;

AOI_FillArray_Real2D(FillLoop,FlowAir,0,2000);
AOI_FillArray_Real2D(FillLoop,PressArr,0,2000);
//ClearArray_2D(FillLoop,FlowO2);           // Clear all raw data

// Main flow

StateMain := 10;
SetupDone := 0;
NewSettings := 1;
StartVent := 0;
StopVent := 0;
VentOn := 0;

// Setup

SetupState := 0;

// GUI

data_changed :=0;

// measCalc

COP(RunModes[0], Mode, 82); // Set Mode to 'VC'
inspStartI := 0;
inspPauseStartI := 0;
expStartI := 0;
expEndI := 0;
totTime := 0; //in s
endI := 0;
PplatBreathFrac := 0.2;

AOI_FillArrayReal50(Fill50, PpeakArr, 0, 50);           //Float, 50 length array with Ppeak from measCalc(), used
for minAvg                                         //to calculate average pressure
AOI_FillArrayReal50(Fill50, PEEPArr, 0, 50);          //Float, 50 length array with positive end expiratory
pressure in cmH2O, from measCalc()                   //pressure
AOI_FillArrayReal50(Fill50, FlowImaxArr, 0, 50);      //Float, 50 length array with max insp flow in l/min, from
measCalc()                                         //measCalc()
AOI_FillArrayReal50(Fill50, PplatArr, 0, 50);         //Float, 50 length array with plateau insp pressure in
cmH2O, from measCalc()                            //cmH2O
AOI_FillArrayReal50(Fill50, RRArr, 0, 50);           //Float, 50 length array with respiratory rate in 1/min,
from measCalc()                                     //from measCalc()
AOI_FillArrayReal50(Fill50, VtiArr, 0, 50);          //Float, 50 length array with insp tidal volume in ml, from
measCalc()                                         //measCalc()
AOI_FillArrayReal50(Fill50, MVArr, 0, 50);           //Float, 50 length array with minute ventilation in l/min,
from measCalc()                                     //from measCalc()
AOI_FillArrayReal50(Fill50, CrsDynArr, 0, 50);        //Float, 50 length array with dynamic respiration system
compliance in ml/cmH2O, from measCalc()             //compliance
AOI_FillArrayReal50(Fill50, CrsQStatArr, 0, 50);      //Float, 50 length array with quasistatic respiratory system
compliance in ml/cmH2O, from measCalc()             //compliance
```

**FirstScan - Structured Text**

test:P\_Sample\_p5:measArray  
Total number of lines in routine: 109

```

AOI_FillArrayReal100(Fill100, RawArr, 0, 50);           //Float, 50 length array with airway resistance in cmH2O/l/ →
s, from measCalc()
AOI_FillArrayReal100(Fill100, BreathDurationArr, 0, 50); //FLOAT, 50 length array with seconds for each breath for →
calculation of minute avgVentStateArr. Calculated from knowing meas at 100Hz for all breath data →

// CalcBreath

StateCalcBreath := 0;          // Reset sequence
valveUpdInterval := 10;        //in ms Frequency 1000/10 = 100Hz
maxFlowAllowed := 100;         //In L/min. Sensor works until 830 ml/s (=49.8 l/min) Note: I will see if we have →
options to work around this as it will severely limit our options for settings →

// See if first time setup is done, the default values else persistent values from controller
if FiO2 = 0 then;
    FiO2 := 4;
else
    FiO2 := FiO2Running;
end_if;

if VT = 0 then;
    VT := 350;
else
    VT := VtRunning;
end_if;

if RRate = 0 then;
    RRate := 15;
else
    RRate := RRateRunning;
end_if;

if IE_I = 0 then;
    IE_I := 1;
else
    IE_I := IE_IRunning;
end_if;

if IE_E = 0 then;
    IE_E := 2;
else
    IE_E := IE_ERunning;
end_if;

if inspPauseTime = 0 then;
    inspPauseTime := 10;
else
    inspPauseTime := inspPauseTimeRunning;
end_if;

if inspRampTime = 0 then;
    inspRampTime := 5;
else
    inspRampTime := inspRampTimeRunning;
end_if;

// Breath detect
_StartPointBefore := 30;
_StartPosition:= 0;

```

```
_endPosition := 0;  
_PointAfterEnd := 30;
```

**General**

Type:  Structured Text      Number of Lines: Original 109  
In Program:  measArray

```
// If 1. scan in PLC
if s:fs then;
    jsr(FirstScan,0);
end_if;

// Jump to subroutine
jsr(measArrays,0);
```

**General**

Type:  
In Program:

 Structured Text  
 measArray

Number of Lines: Original 9

```

// Enables to change scanrate
// Special on a Rockwell PLC were you can change scan rates on the fly
// Settings can be between 50hz or 100hz - alternativ it can be a fixed setting

if Samplerate <> Samplerate_Old then;

    if Samplerate <= 50 then;
        Samplerate := 50;                                // Set samplerate to 50 Hz if less is chosen
    end_if;

    if Samplerate > 100 then;
        Samplerate := 100;                             // Set samplerate to 100 Hz if higher is chosen
    end_if;
    TimeRate := (1000/Samplerate)*1000;           // Calc from Hz to ms
    SSV(Task,P_Sample_p5,Rate,TimeRate);          // Shift sample rate on data collection
    SSV(Task,P_Sample_p8,Rate,TimeRate);          // Shift sample rate on data collection
    Samplerate_Old := Samplerate;                  // Update parameter
    GSV(Task,P_Sample_p5,Rate,NewTimeRate_p5);    // Get the time rate
    GSV(Task,P_Sample_p8,Rate,NewTimeRate_p8);    // Get the time rate

    if TimeRate = NewTimeRate_p5 and TimeRate = NewTimeRate_p8 then;
        SampelRateChangeAllowed := 0;                // Stop change again
        Samplerate_Old := 0;
    end_if;
end_if;

// FIFO fill - load raw analog data into array
if counter >= 0 and counter <= 999 then;          // Test counter values is between array
values
    Flow[counter,0] := (AnalogInputFlow*1)/1000;      // Write Total flow measure into array 1.
dimension
    Flow[counter,0] := 0;                            // Write data read into array 2. dimension
    FlowAir[counter,0] := (AnalogInputFlowAir*75)/1000; // Write flowAir measure into array 1.
dimension
    FlowAir[counter,1] := 0;                         // Write data read into array 2. dimension
    (AnalogInputFlowO2*75)/1000;                      // Write flowO2 measure into array 1.
dimension
    FlowO2[counter,0] := 0;                          // Write data read into array 2. dimension
    (AnalogInputFlowO2*75)/1000;                      // Write pressureArr measure into array 1.
    PressArr[counter,0] := (AnalogInputPressure*6.37322633125)/1000 ; // Write data read into array 2. dimension
    PressArr[counter,1] := 0;                         // Write data read into array 2. dimension

// -----
// ----- C O D E   F R O M   E D B -----
// All my variables start with _edb.
// Filtered values are stored in arrays:
//     _edbFlowArrFiltered
//     _edbFlowO2Filtered
//     _edbPressArrFiltered
// Note that initial setup of filters could be moved in a first scan
//     initialization routine (code between lines:
// **** Setup filters ****
// ...
// **** End of setup filters ****
//
// I'm assuming here that all filter

```

```
// are setup in the same way ( same window size, same type, ... ).  
// In this isn't true the following code must be modified.  
  
if counter = 0 then  
  
    // ***** Setup filters *****  
    // The lines in this section never change. Hence they could be moved  
    // in a first scan initialization.  
  
    // Set the window size of the filter.  
    // The total window size will be: 2 * _edbHalfWindowFilterSize + 1  
    // Maximum allowed value is 15  
    _edbHalfWindowFilterSize := 5; // To be defined  
  
    // Set the type of the filter.  
    // Mean filter (1) is basic a convolution and requires setting the  
    // weights  
    _edbTypeOfFilter := 0; // 0: MEDIAN; 1: MEAN  
    _edbParUseWeights := 0; // Not used with a median filter  
  
    // In case _edbTypeOfFilter = 1 and _edbParUseWeights := 1  
    // weights must be set.  
    // ...  
  
    // After having decided the above parameters, let's setup the  
    // filters:  
    _edbAoiFilterFlowAir.Par_nWindowLeft := _edbHalfWindowFilterSize;  
    _edbAoiFilterFlowAir.Par_nWindowsRight := _edbHalfWindowFilterSize;  
    _edbAoiFilterFlowAir.Par_eType  
    _edbAoiFilterFlowAir.Par_bUseWeights := _edbParUseWeights;  
  
    _edbAoiFilterFlowO2.Par_nWindowLeft := _edbHalfWindowFilterSize;  
    _edbAoiFilterFlowO2.Par_nWindowsRight := _edbHalfWindowFilterSize;  
    _edbAoiFilterFlowO2.Par_eType  
    _edbAoiFilterFlowO2.Par_bUseWeights := _edbParUseWeights;  
    // ***** End of setup filters *****  
  
    // ***** Init filters' cycle *****  
    // The following lines must be execute when an acquisition cycle starts  
    // (they cannot be moved in a first scan initialization routine)  
    _edbAoiFilterFlowAir.In_bInitialize := 1;  
    _edbAoiFilterFlowO2.In_bInitialize := 1;  
    _edbAoiFilterPress.In_bInitialize := 1;  
    // ***** End of init filters' cycle *****  
  
end_if;  
  
// Calculate filtered values  
AOI_FilterSignal( _edbAoiFilterFlowAir, AnalogInputFlowAir, _edbParWeights );  
AOI_FilterSignal( _edbAoiFilterPress, AnalogInputPressure, _edbParWeights );  
AOI_FilterSignal( _edbAoiFilterFlowO2, AnalogInputFlowO2, _edbParWeights );  
  
// Check if initial buffering has finished so we can start storing values.
```

```
// First _edbHalfWindowSize elements of the array won't be calculated.  
// Last _edbHalfWindowSize elements of the array won't be calculated as well.  
_edbCanStore := not _edbAoiFilterFlowAir.Out_Buffering and  
not _edbAoiFilterFlowO2.Out_Buffering and  
not _edbAoiFilterPress.Out_Buffering and  
counter >= 2 * _edbHalfWindowSize; // This last condition should be true when the  
previous conditions are true.  
if _edbCanStore then  
    _edbIdx := counter - _edbHalfWindowSize; // Calculated filtered data is delayed by  
_edbHalfWindowSize acquisitions  
    _edbFlowArrFiltered[_edbIdx,0] := _edbAoiFilterFlowAir.Out_rFilteredSignalDelayed;  
    _edbFlowArrFiltered[_edbIdx,1] := 0;  
    _edbFlowO2Filtered[_edbIdx,0] := _edbAoiFilterPress.Out_rFilteredSignalDelayed;  
    _edbFlowO2Filtered[_edbIdx,1] := 0;  
    _edbPressArrFiltered[_edbIdx,0] := _edbAoiFilterFlowO2.Out_rFilteredSignalDelayed;  
    _edbPressArrFiltered[_edbIdx,1] := 0;  
end_if;  
// ----- E N D O F C O D E F R O M E D B -----  
// -----  
counter := counter +1;                                // next count  
end_if;  
  
if counter >= 999 then;                                // If count is at end of array  
    counter := 0;                                     // Reset counter  
end_if;
```

**General**

Type:  Structured Text      Number of Lines: Original 138  
In Program:  measArray

**General****Configuration**

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

**Monitor****Scan Times(Execution Time)**

Max:	797 us	Last:	37 us
------	--------	-------	-------

```
Pressure := 100.0/16000.0 * (AnalogInputPressure - 4000.0);
//PressHtreshold := ((PressHtresholdUser1/100.0)*16000.0)+4000.0;

// Control that either inspiratory or expiratory is set
if (Insp and Exp) or (not Insp and not Exp) then;
    if inspCount < inspCountLength and expCount < expCountLength then;
        Exp := 1;                                // Select expiatory
        Insp := 0;                               // Deselect inspiratory
    end_if;
    if inspCount < inspCountLength then;
        Insp := 1;                                // Select inspiraty
        Exp := 0;                               // DeSelect expiatory
    end_if;
    if expCount < expCountLength then;
        Insp := 0;                               // Deselect inspiratory
        Exp := 1;                                // Select expiatory
    end_if;
end_if;

if VentOn then;

    if StateValveCtrl = 0 then;

        // Control if Inspiratory or expiatory
        if Insp then;
            StateValveCtrl :=10;                  // Inspiratory
            if _StartPosition <> _StartpositionOld then;
                _StartPosition := counter;
                _StartpositionOld := _StartPosition;
            end_if;
        end_if;

        if Exp then
            StateValveCtrl :=100;                 // Expiatory
        end_if;
    End_if;

    if StateValveCtrl = 10 then;

        // Start of the Inspiratory
        inspCount := inspCount +1;               // Count inspiratory up one in the cycle
        expCount := 0;                          // Reset expiatory counter
input     if Pressure > PressHtresholdUser1 then;           // Analog pressure High warning is high from ▷
            StateValveCtrl := 11;                // Go to high pressure routine
        else
            StateValveCtrl := 20;                // Go to normal routine
        end_if;
    end_if;

    if StateValveCtrl = 11 then;
        // Close valves - high alarm on the Analog input will disappear when pressure is lower than high alarm
        AnalogOutFlowO2 := 0;                  // Close O2 valve
        AnalogOutFlowAir := 0;                 // Close Air Valve
        expOpen :=1;                           // Output - Open expiatory valve
        StateValveCtrl := 30;                 // Next state
    end_if;

    if StateValveCtrl = 20 then;
```

```

// Normal routine - set values to open valve
expOpen := 0;                                // Output - close the expiatory valve
MCSV(GetSlaveValue, O2Up, inSpO2FlowCtrl[inSpCount], O2mAOut, O2mASlope, O2mADer); // Set value to analog output control of O2
AnalogOutFlowO2 := 1.03*O2mAOut;                // Set value to analog output control of Air
MCSV(GetSlaveValue, AirUp, inSpAirFlowCtrl[inSpCount], AirmAOut, AirmASlope, AirmADer); // Set value to analog output control of Air
AnalogOutFlowAir := 1.03*AirmAOut;              // Next state
StateValveCtrl := 30;
end_if;

if StateValveCtrl = 30 then;
    // Control if inspiratory has ended
    if inSpCount >= inSpCountLength then;        // The Inspiratory has ended
        Insp := 0;                                // Deselect inspiratory
        Exp := 1;                                // Select expiatory
        if newInspFlowAvil then;                  // New Inspiratory data is available
            StateValveCtrl := 31;                  // Set to new data for next cyclus is available
        else
            StateValveCtrl := 0;                  // Start over
        end_if;
    else
        StateValveCtrl := 0;                  // Start over
    end_if;
end_if;

if StateValveCtrl = 31 then;
    COP(nextInspO2FlowCtrl[0], inSpO2FlowCtrl[0], 650); // Copy new O2 data to next inspiratory
    COP(nextInspAirFlowCtrl[0], inSpAirFlowCtrl[0], 650); // Copy new Air data to next inspiratory
    := 0;                                              // Reset new Inspiratory data is available
    newInspFlowAvil := 0;                            // Start over
    StateValveCtrl
end_if;

if StateValveCtrl = 100 then;
    // Start of the Expiatory
    expCount := expCount +1;                      // Count expiatory up one in the cycle
    inSpCount := 0;                                // Reset inspiratory counter
    AnalogOutFlowO2 := 0;                          // Close Inspiratory O2 valve
    AnalogOutFlowAir := 0;                         // Close Inspiratory Air valve
    expOpen := expOpenCloseCtrl[expCount].0;       // Set state of expiatory valve
    if _expCountLength then;                      // The Expiratory has ended
        expValveCtrl := 101;
        if _endPosition == _endPositionOld then;
            <= counter;
            _endPositionOld := _endPosition;
        end_if;
    else
        StateValveCtrl := 0;                      // Start over
    end_if;
end_if;

if StateValveCtrl = 101 then;
    Insp := 1;                                // Select inspiratory
    Exp := 0;                                // DeSelect expiatory
    if flowCalcOn then;

```

```
_StartPointBertore := _StartPointBertore;
_StartPosition:= _StartPosition;
_endPosition := _endPosition ;
_PointAtterEnd := _PointAtterEnd;

if (_endPosition - _StartPosition) + _StartPointBertore +_PointAfterEnd < 999 then;

    if _endPosition - _StartPosition > 0 then

        StartSendPosition := _StartPosition - _StartPointBertore;

        if StartSendPosition < 0 then

            _CopSource := 999 + StartSendPosition;
            _CopLength := abs(StartSendPosition+1)*2;
            COP(FlowAir[_CopSource,0],FlowArrMem[0,0],_CopLength);
            COP(FlowO2[_CopSource,0],FlowO2Mem[0,0],_CopLength);
            COP(PressArr[_CopSource,0],PressArrMem[0,0],_CopLength);

            _CopSource1 := 0;
            _CopDes1 := _CopLength + 1;
            _CopLength1 := abs(StartSendPosition)+1;
            COP(FlowAir[_CopSource1,0],FlowArrMem[_CopDes1 ,0],_CopLength1);
            COP(FlowO2[_CopSource1,0],FlowO2Mem[_CopDes1 ,0],_CopLength1);
            COP(PressArr[_CopSource1,0],PressArrMem[_CopDes1 ,0],_CopLength1);

        else

            _CopSource3 := StartSendPosition+1;
            _CopLength3 :=_StartPointBertore*2;
            COP(FlowAir[_CopSource3,0],FlowArrMem[0,0],_CopLength3);
            COP(FlowO2[_CopSource3,0],FlowO2Mem[0,0],_CopLength3);
            COP(PressArr[_CopSource3,0],PressArrMem[0,0],_CopLength3);

        end_if;

        _EndSendPosition := _endPosition + _PointAfterEnd;

        If _EndSendPosition > 999 then

            _CopSource4 := _StartPosition + 1;
            _CopDes4 := _StartPointBertore;
            _CopLength4 := ( 999 - _StartPosition + 1)*2;
            COP(FlowAir[_CopSource4,0],FlowArrMem[_CopDes4,0],_CopLength4);
            COP(FlowO2[_CopSource4,0],FlowO2Mem[_CopDes4,0],_CopLength4);
            COP(PressArr[_CopSource4,0],PressArrMem[_CopDes4,0],_CopLength4);

            _CopSource5 := 0;
            _CopDes5 := _CopDes4 + 999 - _StartPosition;
            _CopLength5 := (_EndSendPosition - 999) *2;
            COP(FlowAir[_CopSource5,0],FlowArrMem[_CopDes5,0],_CopLength5);
            COP(FlowO2[_CopSource5,0],FlowO2Mem[_CopDes5,0],_CopLength5);
            COP(PressArr[_CopSource5,0],PressArrMem[_CopDes5,0],_CopLength5);

        else

            _CopSource6 := _StartPosition;
```

```
_CopDes6 := _StartPointBefore;
_CopLength6 := (_EndSendPosition - _StartPosition) * 2;
COP(FlowAir[_CopSource6,0],FlowArrMem[_CopDes6,0],_CopLength6);
COP(FlowO2[_CopSource6,0],FlowO2Mem[_CopDes6,0],_CopLength6);
COP(PressArr[_CopSource6,0], PressArrMem[_CopDes6,0],_CopLength6);

End_if;

else

    StartSendPosition := _StartPosition - _StartPointBefore;
    _CopSource7 := StartSendPosition;
    _CopLength7 := (999 - StartSendPosition + 1) *2;
    COP(FlowAir[_CopSource7,0],FlowArrMem[0,0],_CopLength7);
    COP(FlowO2[_CopSource7,0],FlowO2Mem[0,0],_CopLength7);
    COP(PressArr[_CopSource7,0], PressArrMem[0,0],_CopLength7);

    _EndSendPosition := _endPosition + _PointAfterEnd;
    _CopSource8 := 0;
    _CopDes8 := (999 - StartSendPosition *2);
    COP(FlowAir[_CopSource8,0],FlowArrMem[_CopDes8,0],_CopLength8);
    COP(FlowO2[_CopSource8,0],FlowO2Mem[_CopDes8,0],_CopLength8);
    COP(PressArr[_CopSource8,0], PressArrMem[_CopDes8,0],_CopLength8);

end_if;
end_if;

_endPositionOld :=0;
_StartpositionOld := 0;
CalcNewData := 1;

else

    CalcNewData := 0;

end_if;

if NewExpOpenCloseAvail then;                                // New Expiatory data is available
    StateValveCtrl := 102;                                    // Next state
else StateValveCtrl := 0;                                    // Start over
end_if;
end_if;

if StateValveCtrl = 102 then;
    COP(nextExpOpenCloseCtrl[0], expOpenCloseCtrl[0], 1000); // Copy new close/open data to next expiatory

    expCountLength := nextExpCountLength;                    // Copy new lenght of expiatory
    NewExpOpenCloseAvail := 0;                                // Reset Expiatory data is available
    StateValveCtrl := 0;                                    // Start over
end_if;
end_if;
```

**General**

Type: Structured Text      Number of Lines: Original 229  
In Program: ISR\_ValveCtrl

**General****Configuration**

Type:	Periodic	Watchdog:	500.000 ms
Period:	10.000 ms	Disable automatic output processing to reduce task overhead:	No
Priority:	8	Inhibit task:	No

**Program Schedule****Scheduled**

Setup	calcBreath	initVent	BreathDetect
measCalc	MainFlow		

**Unscheduled****Monitor****Scan Times(Elapsed Time)**

Max:	4.712000 ms	Last:	1.271000 ms
------	-------------	-------	-------------

**Interval Times(Elapsed Times Between Triggers)**

Max:	10.851000 ms	Min:	9.203000 ms
------	--------------	------	-------------

Task overlap count:

0

**General****Configuration**

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

**Monitor****Scan Times(Execution Time)**

Max:	37 us	Last:	25 us
------	-------	-------	-------

```
// Init program of the GUI (HMI panel)
// Write initial code here
```

**General**

Type:  Structured Text      Number of Lines: Original 6  
In Program:  Setup

```
// Init program of the ISR_RS232 AI controller  
// Write initial coder here
```

**General**

Type:  Structured Text      Number of Lines: Original 5  
In Program:  Setup

```
// Init program of the Init variables for the program  
// Write initial coder here
```

**General**

Type:  Structured Text      Number of Lines: Original 3  
In Program:  Setup

```
case SetupState of
  10: JSR(InitVariables,0);           // Jump to Init Variables
      SetupState := 20;               // Set next state

  20: JSR(InitGUI,0);                // Jump to Init GUI
      SetupState := 30;               // Set next state

  30: JSR(InitISR_RS232);          // Jump to init ISR_RS232
      SetupDone :=1;                 // Set Setup Done
      SetupState := 99;               // Sequence Done

end_case;
```

**General**

Type:  Structured Text      Number of Lines: Original 13  
In Program:  Setup

**General****Configuration**

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

**Monitor****Scan Times(Execution Time)**

Max:	2606 us	Last:	8 us
------	---------	-------	------

```
//Saves o2 and air inspiratory flows in l/min as 100 Hz arrays for ISR_valveCtrl() to
//read and convert to a current securing these defined flows delivered to
//the patient.

if StateCalcBreath = 10 then;

//Calculate durations
inspPauseTimePercent := inspPauseTime/100; //Percentage of whole inspIE_In where there is a pause (in the end)
inspRampTimePercent := inspRampTime/100; //Percentage of active inspIE_In (when there is flow)
BreathDurationInMS := (1/RRate) * 60000;
inspTimeinMS := IE_I * (BreathDurationInMS/(IE_I+IE_E));
totalCount := BreathDurationInMS/valveUpdInterval;
inspCounts := inspTimeinMS/valveUpdInterval;
AOI_Floor(Floor,inspCounts,inspCountLengthNew);
expTimeInMS := BreathDurationInMS - inspTimeinMS;
expCountlengthCalc := totalCount-inspCountLengthNew;
AOI_Floor(Floor, expCountlengthCalc, expCountLengthNew);

//Calculate part of inspIE_In where valves are open
inspPauseTimeInMS := inspPauseTimePercent * BreathDurationInMS;
InspiratoryRampTimeInMS := inspRampTimePercent * BreathDurationInMS;
ConstantFlowTimeInMS := inspTimeinMS - inspPauseTimeInMS - InspiratoryRampTimeInMS;

//Calculate flow in L/min required to reach Vt in active inspIE_In
//(assuming a linear ramp)
ConstantInspiratoryFlow := 1/(ConstantFlowTimeInMS+0.5*InspiratoryRampTimeInMS)*Vt*60;

//This needs to be warned to the user - We will need a handle in HMI to not
//allow such settings if we are stuck with this limit on max flow
if ConstantInspiratoryFlow > maxFlowAllowed then;
    COP(Messages_to_user[2], message, 82);
end_if;

//Calculate split between O2 and Air
FractionAir := ((1-FiO2/100)/0.79);
FractionO2 := 1 - FractionAir;

//Calculate Ramp samples for O2 and Air
O2ConstantFlow := FractionO2 * ConstantInspiratoryFlow; //L/min
AirConstantFlow := FractionAir * ConstantInspiratoryFlow; //L/min
RampNrCalc := InspiratoryRampTimeInMS/valveUpdInterval;
AOI_Floor(Floor, RampNrCalc,RampNrOfAdjustments);
O2RampFlowIncreasePerAdjust := O2ConstantFlow/RampNrOfAdjustments;
AirRampFlowIncreasePerAdjust := AirConstantFlow/RampNrOfAdjustments;

ConstansFlowNrCalc := ConstantFlowTimeInMS/valveUpdInterval;
AOI_Floor(Floor, ConstansFlowNrCalc,ConstantFlowNrOfAdjustments );
    inspPauseTimeInMS/valveUpdInterval;

inspPauseCalc :=
AOI_Floor(Floor, inspPauseCalc,inspPauseNrOfAdjustments );
//Build control arrays
AOI_FillArrayReal(Fill, inspO2FlowCtrlNew, 0, 650);
AOI_FillArrayReal(Fill, inspAirFlowCtrlNew, 0, 650);
AOI_FillArrayInt(Fillint, expOpenCloseCtrlNew, 0, 1000);

//inspiratory
for i :=0 to inspCountLengthNew by 1 do;
```

```
if i <= RampNrOfAdjustments then;
    RampTime[i] := valveUpdInterval*(i-1);
    inspO2FlowCtrlNew[i] := i * O2RampFlowIncreasePerAdjust;
    inspAirFlowCtrlNew[i] := i * AirRampFlowIncreasePerAdjust;
elseif i > RampNrOfAdjustments and i < RampNrOfAdjustments+ConstantFlowNrOfAdjustments then;
    inspO2FlowCtrlNew[i] := O2ConstantFlow;
    inspAirFlowCtrlNew[i] := AirConstantFlow;
elseif i > (RampNrOfAdjustments+ConstantFlowNrOfAdjustments) and i <= (RampNrOfAdjustments
+ConstantFlowNrOfAdjustments+inspPauseNrOfAdjustments) then;
    inspO2FlowCtrlNew[i] := 0;
    inspAirFlowCtrlNew[i] := 0;
end_if;
end_for;

//expiratory
AOI_FillArrayInt(Fillint, expOpenCloseCtrlNew, 1, expCountLengthNew);

// Copy calculated data for Inspiratoty to next breath cycles
COP(inspO2FlowCtrlNew[0], nextInspO2FlowCtrl[0], 650);
COP(inspAirFlowCtrlNew[0], nextInspAirFlowCtrl[0], 650);
nextInspCountLength := inspCountLengthNew;

// Copy calculated data for Expiratory to next breath cycles
COP(expOpenCloseCtrlNew[0], nextExpOpenCloseCtrl[0], 1000);
nextExpCountLength := expCountLengthNew;

NewSettings := 0;
newInspFlowAvil := 1;
NewExpOpenCloseAvail:= 1;
StateCalcBreath := 0;

end_if;
```

**General**

Type:  Structured Text      Number of Lines: Original 92  
In Program:  calcBreath

**General****Configuration**

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

**Monitor****Scan Times(Execution Time)**

Max:	195 us	Last:	10 us
------	--------	-------	-------

```
// Init of Valve Control
// Routine only runs once

// Fill BreathDurationArr with -1 only on first breath, for messuring average values
AOI_FillArrayReal50(Fill50,BreathDurationArr, -1, 50);
flowCalcOn := 0;           // No calculation of measurement - set to 1 if measurement ***** TEMP Soulution *****
```

**General**

Type: Structured Text      Number of Lines: Original 8  
In Program: initVent

```
// Init of Pressure Watch  
// Routine only runs once
```

**General**

Type:  
In Program:

 Structured Text  
 initVent

Number of Lines: Original 2

```
case StateInitVent of

    10:                                // Prepare for transfer of data from
        calcBreath routine

        COP(nextInspO2FlowCtrl[0], inspO2FlowCtrl[0], 650);
        COP(nextInspAirFlowCtrl[0], inspAirFlowCtrl[0], 650);
        COP(nextExpOpenCloseCtrl[0], expOpenCloseCtrl[0], 1000);
        inspCountLength := nextInspCountLength;
calculated to running data
        expCountLength := nextExpCountLength;
calculated to running data

        if inspCountLength >20 and expCountLength >20 then;
            StateInitVent := 20;
        else
            StartVent :=0;
            VentOn := 0;
            COP(Messages_to_user[1], message, 82);
not ok
        end_if;

    20:
        JSR(initISR_ValveCtrl);
runs once
        JSR(Init_ISR_PressWatchID);
runs once

        VentOn := 1;
        StartVent := 0;
        StateInitVent := 99;

end_case;
```

// Copy calculated data to running data  
// Copy calculated data to running data  
// Copy calculated data to running data  
// Copy lenght of inspiratory array  
  
// Copy lenght of expiratory array  
  
// Check values is correct  
// Set next stage  
  
// Stop ventilation  
// Set do not run to ISR\_ValveCtrl  
// Write to user that lenght insp or exp is

// Jump to init for ISR\_ValveCtrl - only  
// Jump to init for ISR\_PressWatchID - only

// Set start to ISR\_ValveCtrl  
// Reset start bit  
// Set sequence finished state

**General**

Type:  Structured Text      Number of Lines: Original 28  
In Program:  initVent

## General

### Configuration

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

## Monitor

### Scan Times(Execution Time)

Max:	37 us	Last:	7 us
------	-------	-------	------



**General**

Type:  
In Program:

 Structured Text  
 BreathDetect

Number of Lines: Original 1

**General****Configuration**

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

**Monitor****Scan Times(Execution Time)**

Max:	36 us	Last:	18 us
------	-------	-------	-------

```
// from inflow and press alone (may require filtering such as mean/median!)
```

```
// Only check indices with FlowArr[0,x]=0, meaning not analysed yet
```

```
if CalcNewData then;
```

```
    inspStartI := 0;
```

```
    expStartI := 0;
```

```
    expEndI := 0;
```

```
    for _i := 1 to 999 by 1 do
```

```
        if FlowArrMem[_i,1] = 0 then; // Check if data is read already
```

```
            if inspStartI = 0 then; // If Inspiratory start point ↗
```

```
                not found // Start of insp is at shift from flow of zero to flow of >0.
```

```
                if FlowArrMem[_i-1,0] = 0 and FlowArrMem[_i,0] > 0 then; // If start of climbing flow, ↗
```

```
were last point was zero and this point has a value
```

```
                    inspStartI := _i ; // Set Inspiratory start point ↗
```

```
                    end_if;
```

```
                    end_if;
```

```
                if inspStartI > 0 and expStartI = 0 then; // If Inspiratory start point ↗
```

```
found but nor the Expiratory point
```

```
                    // Start of exp is at shift from to flow of >0 to flow of zero.
```

```
                    if FlowArrMem[_i-1,0] > 0 and FlowArrMem[_i,0] = 0 then; // If end of climbing flow, and ↗
```

```
start of Pplat + Pplat time
```

```
                        if Mode = RunModes[0] then; // If mode is VC
```

```
                            inspPauseStartI := _i; // Set inspiratory pause start ↗
```

```
point
```

```
                        expStartI := _i + ABS(tipause/(TimeRate/1000)); // Set Expiratory start point ↗
```

```
                        else // Set Expiratory start point ↗
```

```
                            expStartI := _i;
```

```
                            end_if;
```

```
                        end_if;
```

```
                    end_if;
```

```
                if inspStartI > 0 and expStartI > 0 and expEndI = 0 then; // If Inspiratory start point ↗
```

```
found and the Expiratory start point, but not Expiratory end point
```

```
                    // End of exp is at shift from to flow of zero to flow of >0.
```

```
                    if FlowArrMem[_i-1,0] = 0 and FlowArrMem[_i,0] > 0 then; // If flow is zero and detection ↗
```

```
of new breath
```

```
                        expEndI := _i - 1; // Set Expiratory end point ↗
```

```
                        end_if;
```

```
                    end_if;
```

```
    end_for;
```

```
// Check of illegal calculation that would bring the controller in Major fault
```

```
if inspStartI < 1 then;
```

```
    inspStartI := 1 ;
```

```
end_if;
```

```

if expStartI <1 then;
    expStartI := inspStartI +1;
end_if;

if expEndI < 1 then;
    expEndI := expStartI +1;
end_if;

// Set all calculations vaerables for breath

// Breathduration
BreathDuration := ((expEndI-inspStartI+1)/Samplerate);
RR := (60/BreathDuration);

// in area with inflow=0 before insp start --> PEEP
PEEP := PressArrMem[inspStartI-1,0];

// in area with inflow=0 check for increase in press (Exp pause, exp flow=0) --> PEEP+PEEPI
// NOTE as exp pause cannot be detected this is not possible at the moment,
// but should be done here
//PEEP := (PEEP+PEEPI)-PEEP;

// Peak = Max pressure during insp
// NOTE: will require filtering, e.g. moving average

AOI_Max(PpeakMaxValue,PressArrMem, 0, inspStartI, expStartI);
Ppeak := PpeakMaxValue.Out_Max;

// Flowi,max = max flow during insp
// NOTE: will require filtering e.g. moving average, should be around same
// sample range as Ppeak
// in VC this is the constant delived level of insp Flow

AOI_Max(FlowImaxValue, FlowArrMem, 0,inspStartI, expStartI);
FlowImax := FlowImaxValue.Out_Max;

// Pplat = pressure at period immediately following inFlow=0 where
// pressure equilibrates
// only available if insp pause was performed
// Calculated as mean of last fifth of the pause period

        then;
        (inspPauseTimeInMS/1000) * Samplerate*PplatBreathFrac;
if inspPauseStartT > 0
    AOI_Round(Round,StartPointValue, StartPointMean);
    StartPointValue := expStartI-StartPointMean-1;
    EndPoint := expStartI-1;
AOI_Mean2D(PplatMean, PressArrMem, 0, StartPoint, EndPoint);
    Pplat := PplatMean.Out_Mean;
else
    Pplat := 0;
end_if;

//Vti = integrate inFlow(insp start to exp start)

eod_for;= inspStartI to expStartI by 1 do
    Vti := Vti + FlowArrMem[_J,0] * (1/(Samplerate*60));
//MV = RR*Vti

```

```
MV := RR*Vti;

//Raw = (Ppeak-Pplat)/Flowi,max    %NOTE: in cmH2O/l/s (different time
//unit than Flowi,max at l/min
//only available if insp pause was performed so that we have Pplat
//assumes FlowIMax is at PIP

if inspPauseStartI > 0 then;
    Raw := (Ppeak-Pplat)/(FlowIMax/60); //FlowIMax/60 to go from l/min to l/s as Raw is in units of l/s
else
    Raw := 0;
end_if;

//Crs_dyn = Vt/(Ppeak-PEEP)
CrsDyn := Vti/(Ppeak-PEEP);

//Crs_qstat = Vt/(Pplat-PEEP)
if inspPauseStartI > 0 then;
    CrsQStat := Vti/(Pplat-PEEP);
else
    CrsQStat := 0;
end_if;

//set second dimension status marker in FlowArr and PressArr to 1 from
//inspStartI to expEndI to indicate that this breath has been analysed
for _n := inspStartI to expEndI+1 by 1 do

    //FlowArrMem[_n,1] := 1;
    //PressArrMem[_n,1] := 1;

end_for;

for _T := 49 to 1 by -1 do

    PpeakArr[_T] := PpeakArr[_T-1];
    PEEPArr[_T] := PEEPArr[_T-1];
    FlowIMaxArr[_T] := FlowIMaxArr[_T-1];
    PplatArr[_T] := PplatArr[_T-1];
    RRArr[_T] := RRArr[_T-1];
    VtiArr[_T] := VtiArr[_T-1];
    MVArr[_T] := MVArr[_T-1];
    CrsDynArr[_T] := CrsDynArr[_T-1];
    CrsQStatArr[_T] := CrsQStatArr[_T-1];
    RawArr[_T] := RawArr[_T-1];
    BreathDurationArr[_T] := BreathDurationArr[_T-1];

end_for;

    := Ppeak;
PpeakArr[0]:= PEEP;
FlowIMaxArr[0] := FlowIMax;
    := Pplat;
RRArr[0] RR;
VtiArr[0] := Vti;
MVArr[0] := MV;
CrsDynArr[0] := CrsDyn;
```

```

CrsQStatArr[0] := CrsQStat;
RawArr[0] := Raw;
BreathDurationArr[0] := BreathDuration;

//in case there is not a minute of measurements available yet, then calculate
//from the available entries in array
//This could be signaled by setting all values in array as -1 until they
//have been used, this is done by initVent
//I assume here that if an entry is -1 then it is not available

for _z := 1 to 49 by 1 do
  if BreathDurationArr[_z] = -1 and _z = 1 then;
    //disp('error, no breath by breath measurements available');
    _z :=50;                                // Break the for loop
  end_if;
  if _z < 50 then;
    if BreathDurationArr[_z] = -1 then;
      endI := i-1;
      _z :=50;
    end_if;
  end_if;
end_if;

if _z < 49 then;
  totTime := totTime+BreathDurationArr[_z];
end_if;
if totTime >= 60 then;
  endI := _z;
  _z :=50;
end_if;
end_for;

// Calculate the average from last minute

AOI_Mean50(Mean50, PpeakArr, 1, endI);           //1 min avg of peak insp pressure in cmH20
avgPpeak :=Mean50.Out_Mean;                         // Output of mean
AOI_Mean50(Mean50, PEEPArr, 1, endI);             //1 min avg of positive end expiratory pressure in cmH20
avgPEEP :=Mean50.Out_Mean;                          // Output of mean
AOI_Mean50(Mean50, FlowImaxArr, 1, endI);        //1 min avg of max insp flow in l/min
avgFlowImax :=Mean50.Out_Mean;                      // Output of mean
AOI_Mean50(Mean50, PplatArr, 1, endI);            //1 min avg of plateau insp pressure in cmH20
avgPplat :=Mean50.Out_Mean;                         // Output of mean
AOI_Mean50(Mean50, RRArr, 1, endI);               //1 min avg of respiratory rate in 1/min
avgRR :=Mean50.Out_Mean;                            // Output of mean
AOI_Mean50(Mean50, VtiArr, 1, endI);              //1 min avg of insp tidal volume in ml
avgVti :=Mean50.Out_Mean;                           // Output of mean
AOI_Mean50(Mean50, MVArr, 1, endI);               //1 min avg of minute ventilation in l/min
avgMV :=Mean50.Out_Mean;                            // Output of mean
AOI_Mean50(Mean50, CrsDynArr, 1, endI);           //1 min avg of dynamic respiratory system compliance in ml/cmH20
avgCrsDyn :=Mean50.Out_Mean;                        // Output of mean
AOI_Mean50(Mean50, CrsQStatArr, 1, endI);          //1 min avg of quasistatic respiratory system compliance in ml/
cmH20                                            // Output of mean

avgCrsQStat :=Mean50.Out_Mean;                      // Output of mean
AOI_Mean50(Mean50, RawArr, 1, endI);               //1 min avg of airway resistance in cmH20/l/s
avgRaw :=Mean50.Out_Mean;                           // Output of mean

CalcNewData :=0;
AOI_FillArray_Real2D(FillLoop, FlowArrMem, 0, 2000); // Clear all raw data

```

```
AOI_FillArray_Real2D(FillLoop,FlowO2Mem ,0,2000);           // Clear all raw data
AOI_FillArray_Real2D(FillLoop,PressArrMem,0,2000);           // Clear all raw data
end_if;
```

**General**

Type:  Structured Text      Number of Lines: Original 232  
In Program:  measCalc

**General****Configuration**

Main:	 Main	Inhibit program:	No
Fault:	None	Synchronous redundancy data after execution :	Yes

**Monitor****Scan Times(Execution Time)**

Max:	937 us	Last:	30 us
------	--------	-------	-------

```
// Set MainFlow to start over
if StateMain = 0 then;
    StateMain := 10;
end_if;

IE_E := IE_counter +1;

// see if new data is valible from GUI
CsumData := FiO2+Vt + RRate + IE_I + IE_E + inspPauseTime + inspRampTime;

If CsumData <> CsumDataOld then;
    CsumDataOld := CsumData;
    NewSettings :=1;
    FiO2Running := FiO2;
    VtRunning := Vt;
    RRateRunning := RRate;
    IE_IRunning := IE_I;
    IE_ERunning := IE_E;
    inspPauseTimeRunning := inspPauseTime;
    inspRampTimeRunning := inspRampTime;
end_if;

// MainFlow sequence
case StateMain of

    // Run Setup
    10: if SetupDone then;                                // If setup is done
        StateMain := 20;                                 // Set next stage
        SetupState := 0;
    else
        if SetupState = 0 then;
            SetupState := 10;                            // Start the setup program
        end_if;
    end_if;

    // Run CalcBreath
    20: if NewSettings then;                            // If changes has happend on the GUI
        if StateCalcBreath = 0 then;                    // Start the CalcBreath setup
            StateCalcBreath := 10;
        end_if;
    else
        StateMain := 30;                               // Set next stage
    end_if;

    // Stop Ventilation
    30: if StopVent then;                            // Stop ventilation
        if Exp then
            if expCount
                expiratory ends > (expCountLength - 20) then;
                    VentOn :=0;                           // Stop ventilation 20*scanrate befor
                                                    // Stop ventilation routine
                    AOI_FillArray_Real2D(FillLoop,FlowAir,0,2000); // Clear all raw data
                    AOI_FillArray_Real2D(FillLoop,FlowO2,0,2000); // Clear all raw data
                    AOI_FillArray_Real2D(FillLoop,PressArr,0,2000); // Clear all raw data
    end_if;
```

```
AOI_FillArray_Real2D(FillLoop,FlowArrMem,0,2000);           // Clear all raw data
AOI_FillArray_Real2D(FillLoop,FlowO2Mem ,0,2000);           // Clear all raw data
AOI_FillArray_Real2D(FillLoop,PressArrMem,0,2000);           // Clear all raw data
AOI_FillArray_Real2D(FillLoop,_edbFlowArrFiltered,0,2000);   // Clear all raw data
AOI_FillArray_Real2D(FillLoop,_edbFlowO2Filtered,0,2000);   // Clear all raw data
AOI_FillArray_Real2D(FillLoop,_edbPressArrFiltered,0,2000); // Clear all raw data
StopVent := 0;                                              // Reset the stop command
Insp := 1;                                                   // Select inspiratory next start first
Exp := 0;                                                   // Deselect expiratory next start
expCount := 0;                                              // Reset the expiratory count
StateMain := 0;                                              // Reset sequence and start over

    end_if;
end_if;
else
    StateMain := 40;                                         // Set next stage
end_if;

// Run initVent
40: if StartVent then;                                     // If ventilator is started

    if StateInitVent = 0 then;
        StateInitVent := 10;                                 // Start the initVent program
    end_if;

else
    StateMain := 0;
    StateInitVent :=0;

end_if;

end_case;
```

**General**

Type:  
In Program:

 Structured Text  
 MainFlow

Number of Lines: Original 91


**Signature Listing**

**Available Languages**

Relay Ladder



Function Block



Structured Text

```
AOI_ARRAYSHIFT(Inp_Array, _Array);
```

**Parameters**

Required	Name	Data Type	Usage	Description
X	AOI_ARRAYSHIFT	AOI_ARRAYSHIFT	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	Inp_Array	REAL	InOut	
X	_Array	REAL	InOut	

**Extended Description****Execution****Condition      Description**

EnableIn is true

**Revision v1.0 Notes**

Name	Default	Data Type	Scope
Inp_Array		REAL[1000,2]	AOI_ARRAYSHIFT
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Inp_Array - AOI_ARRAYSHIFT/Logic - *0(COP), 0(COP)</i>			
_Array		REAL[1000,2]	AOI_ARRAYSHIFT
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>_Array - AOI_ARRAYSHIFT/Logic - *0(COP), 0(COP)</i>			

Name	Default	Data Type	Scope
<b>No Tags Exist</b>			

**AOI\_ARRAYSHIFT Instruction Definition - Logic Routine**


test:Add-On Instructions:AOI\_ARRAYSHIFT:Logic

Total number of rungs in routine: 1

&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_ARRAYSHIFT


0




(End)

## Available Languages

Relay Ladder



Function Block



Structured Text

```
AOI_Ceil(In_X, Out_Ret);
```

## Parameters

Required	Name	Data Type	Usage	Description
X	AOI_Ceil	AOI_Ceil	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_X	REAL	Input	
X	Out_Ret	REAL	InOut	

## Extended Description

## Execution

Condition	Description
EnableIn is true	

## Revision v1.1 Notes


Name	Default	Data Type	Scope
<b>In_X</b>	0.0	REAL	AOI_Ceil
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_X - AOI_Ceil/Logic - #1, #3</i>			
<b>Out_Ret</b>	??	REAL	AOI_Ceil
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Out_Ret - AOI_Ceil/Logic - *#6</i>			

Name	Default	Data Type	Scope
<u><b>_aoiTrunc</b></u>		AOI_Trunc	AOI_Ceil
Usage:	Local Tag		
External Access:	None		
$_aoiTrunc - AOI\_Ceil/Logic - *#1$			
<u><b>_aoiTrunc.EnableIn</b></u>	1	BOOL	
Enable Input - System Defined Parameter			
<u><b>_aoiTrunc.EnableOut</b></u>	0	BOOL	
Enable Output - System Defined Parameter			
<u><b>_n</b></u>	0	DINT	AOI_Ceil
Usage:	Local Tag		
External Access:	None		
$_n - AOI\_Ceil/Logic - #3, #4, #6, *#2, *#4$			
<u><b>_trunc</b></u>	0.0	REAL	AOI_Ceil
Usage:	Local Tag		
External Access:	None		
$_trunc - AOI\_Ceil/Logic - #2, *#1$			


```
AOI_Trunc( _aoiTrunc, In_X, _trunc );
_n := _trunc;
if _n < In_X then
  _n := _n + 1;
end_if;
Out_Ret := _n;
```

**Available Languages**

Relay Ladder



Function Block



Structured Text

AOI\_FillArrayInt(Array, Source, Len);

**Parameters**

Required	Name	Data Type	Usage	Description
X	AOI_FillArrayInt	AOI_FillArrayInt	In/Out	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	Array	INT	In/Out	
X	Source	DINT	Input	
X	Len	DINT	Input	

**Extended Description****Execution****Condition      Description**

EnableIn is true

**Revision v1.0 Notes**

Name	Default	Data Type	Scope
<b>Array</b>		INT[1000]	AOI_FillArrayInt
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Array - AOI_FillArrayInt/Logic - *0(FLL)</i>			
<b>Len</b>	0	DINT	AOI_FillArrayInt
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Len - AOI_FillArrayInt/Logic - 0(FLL)</i>			
<b>Source</b>	0	DINT	AOI_FillArrayInt
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Source - AOI_FillArrayInt/Logic - 0(FLL)</i>			

Name	Default	Data Type	Scope
<b>No Tags Exist</b>			

**AOI\_FillArrayInt Instruction Definition - Logic Routine**

test:Add-On Instructions:AOI\_FillArrayInt:Logic

Total number of rungs in routine: 1

&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArrayInt


0

FLL	Source	Source
Dest	Array[0]	Len

(End)

## Available Languages

Relay Ladder



Function Block



Structured Text

AOI\_FillArrayReal(Array, Source, Len);

## Parameters

Required	Name	Data Type	Usage Description
X	AOI_FillArrayReal	AOI_FillArrayReal	InOut
	EnableIn	BOOL	Input
	EnableOut	BOOL	Output
X	Array	REAL	InOut
X	Source	DINT	Input
X	Len	DINT	Input

## Extended Description

## Execution

Condition	Description
EnableIn is true	

## Revision v1.0 Notes

**AOI\_FillArrayReal Instruction Definition - Parameter Listing**

test:Add-On Instructions:AOI\_FillArrayReal

Data Type Size: 12 byte (s)

&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArrayReal

Name	Default	Data Type	Scope
<b>Array</b>		REAL[650]	AOI_FillArrayReal
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Array - AOI_FillArrayReal/Logic - *0(FLL)</i>			
<b>Len</b>	0	DINT	AOI_FillArrayReal
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Len - AOI_FillArrayReal/Logic - 0(FLL)</i>			
<b>Source</b>	0	DINT	AOI_FillArrayReal
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Source - AOI_FillArrayReal/Logic - 0(FLL)</i>			

Name	Default	Data Type	Scope
<b>No Tags Exist</b>			


**AOI\_FillArrayReal Instruction Definition - Logic Routine**

test:Add-On Instructions:AOI\_FillArrayReal:Logic

Total number of rungs in routine: 1


&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArrayReal




**Available Languages**

Relay Ladder



Function Block



Structured Text

```
AOI_FillArrayReal50(Array, Source, Len);
```

**Parameters**

Required	Name	Data Type	Usage	Description
X	AOI_FillArrayReal50	AOI_FillArrayReal50	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	Array	REAL	InOut	
X	Source	DINT	Input	
X	Len	DINT	Input	

**Extended Description****Execution****Condition      Description**

EnableIn is true

**Revision v1.0 Notes**

**AOI\_FillArrayReal50 Instruction Definition - Parameter Listing**

test:Add-On Instructions:AOI\_FillArrayReal50

Data Type Size: 12 byte (s)

&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArrayReal50

Name	Default	Data Type	Scope
<b>Array</b>		REAL[50]	AOI_FillArrayReal50
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Array - AOI_FillArrayReal50/Logic - *0(FLL)</i>			
<b>Len</b>	0	DINT	AOI_FillArrayReal50
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Len - AOI_FillArrayReal50/Logic - 0(FLL)</i>			
<b>Source</b>	0	DINT	AOI_FillArrayReal50
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Source - AOI_FillArrayReal50/Logic - 0(FLL)</i>			

Name	Default	Data Type	Scope
<b>No Tags Exist</b>			


**AOI\_FillArrayReal50 Instruction Definition - Logic Routine**

test:Add-On Instructions:AOI\_FillArrayReal50:Logic

Total number of rungs in routine: 1

&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArrayReal50




## Available Languages

### Relay Ladder



### Function Block



### Structured Text

```
AOI_FillArray_Real2D(Array, Source, Len);
```

## Parameters

Required	Name	Data Type	Usage	Description
X	AOI_FillArray_Real2D	AOI_FillArray_Real2D	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	Array	REAL	InOut	
X	Source	DINT	Input	
X	Len	DINT	Input	

## Extended Description

## Execution

### Condition      Description

EnableIn is true

## Revision v1.0 Notes

**AOI\_FillArray\_Real2D Instruction Definition - Parameter Listing**

test:Add-On Instructions:AOI\_FillArray\_Real2D

Data Type Size: 12 byte (s)

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArray\_Real2D

&lt;definition&gt;

Name	Default	Data Type	Scope
<b>Array</b>		REAL[1000,2]	AOI_FillArray_Real2D
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Array - AOI_FillArray_Real2D/Logic - *0(FLL)</i>			
<b>Len</b>	0	DINT	AOI_FillArray_Real2D
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Len - AOI_FillArray_Real2D/Logic - 0(FLL)</i>			
<b>Source</b>	0	DINT	AOI_FillArray_Real2D
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>Source - AOI_FillArray_Real2D/Logic - 0(FLL)</i>			

Name	Default	Data Type	Scope
<b>No Tags Exist</b>			


**AOI\_FillArray\_Real2D Instruction Definition - Logic Routine**

test:Add-On Instructions:AOI\_FillArray\_Real2D:Logic

Total number of rungs in routine: 1


&lt;definition&gt;

LogixProgram\_L27\_Ver4.ACD Data Context: AOI\_FillArray\_Real2D




## Available Languages

## Relay Ladder



## Function Block



## Structured Text

```
AOI_FilterSignal(In_rRawSignal, Par_rWeights);
```

## Parameters

Required	Name	Data Type	Usage	Description
X	AOI_FilterSignal	AOI_FilterSignal	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_rRawSignal	REAL	Input	
	In_bInitialize	BOOL	Input	
	Par_eType	DINT	Input	0: MEDIAN; 1: MEAN
	Par_nWindowLeft	DINT	Input	
	Par_nWindowsRight	DINT	Input	
	Par_bUseWeights	BOOL	Input	Used for mean filtering
	Par_rWeights	REAL	InOut	Used for mean filtering
	Out_Buffering	BOOL	Output	
	Out_rFilteredSignalDelayed	REAL	Output	
	ER	BOOL	Output	
	Err	DINT	Output	

## Extended Description

## Execution

### Condition      Description

EnableIn is true

Prescan

## Revision v1.1 Notes

**AOI\_FilterSignal Instruction Definition - Parameter Listing**

test:Add-On Instructions:AOI\_FilterSignal

Data Type Size: 468 byte (s)

Data Context: AOI\_FilterSignal &lt;definition&gt;

Name	Default	Data Type	Scope
<b>ER</b>	0	BOOL	AOI_FilterSignal
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read Only		
<i>ER - AOI_FilterSignal/Logic - #1, *#19, *#7</i>			
<b>Err</b>	0	DINT	AOI_FilterSignal
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Err - AOI_FilterSignal/Logic - #18, *#14, *#15, *#16, *#17, *#8</i>			
<b>In_bInitialize</b>	0	BOOL	AOI_FilterSignal
Usage:	Input Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>In_bInitialize - AOI_FilterSignal/Logic - #1, #6, *#46</i>			
<b>In_rRawSignal</b>	0.0	REAL	AOI_FilterSignal
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_rRawSignal - AOI_FilterSignal/Logic - #49</i>			
<b>Out_Buffering</b>	0	BOOL	AOI_FilterSignal
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Out_Buffering - AOI_FilterSignal/Logic - #53, *#51</i>			
<i>Out_Buffering - AOI_FilterSignal/Prescan - *#14</i>			
<b>Out_rFilteredSignalDelayed</b>	0.0	REAL	AOI_FilterSignal
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read Only		
<i>Out_rFilteredSignalDelayed - AOI_FilterSignal/Logic - *#75, *#77, *#90</i>			
<b>Par_bUseWeights</b>	0	BOOL	AOI_FilterSignal
Used for mean filtering			
Usage:	Input Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Par_bUseWeights - AOI_FilterSignal/Logic - #25</i>			
<b>Par_eType</b>	0	DINT	AOI_FilterSignal
0: MEDIAN; 1: MEAN			
Usage:	Input Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Par_eType - AOI_FilterSignal/Logic - #17, #17, #31</i>			
<b>Par_nWindowLeft</b>	0	DINT	AOI_FilterSignal
Usage:	Input Parameter		
Required:	No		

**Par\_nWindowLeft (Continued)**

Visible: Yes

External Access: Read/Write

*Par\_nWindowLeft - AOI\_FilterSignal/Logic - #13, #14, #26***Par\_nWindowsRight**

0

DINT

AOI\_FilterSignal

Usage: Input Parameter

Required: No

Visible: Yes

External Access: Read/Write

*Par\_nWindowsRight - AOI\_FilterSignal/Logic - #13, #15, #27***Par\_rWeights**

REAL[1]

AOI\_FilterSignal

Used for mean filtering

Usage: InOut Parameter

Required: Yes

Visible: Yes

Constant: No

*Par\_rWeights - AOI\_FilterSignal/Logic - #32**Par\_rWeights - AOI\_FilterSignal/Prescan - #7*

Name	Default	Data Type	Scope
_bFirstScan	0	BOOL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_bFirstScan - AOI_FilterSignal/Logic - #1, #6, *#45			
_bFirstScan - AOI_FilterSignal/Prescan - *#11			
_bOdd	0	BOOL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_bOdd - AOI_FilterSignal/Logic - #74, *#30			
_bW	0	BOOL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_bW - AOI_FilterSignal/Logic - #37, #87, *#25			
_eType	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_eType - AOI_FilterSignal/Logic - #55, *#31			
_i	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_i - AOI_FilterSignal/Logic - #40, #87, #88, *#39, *#86			
_idx	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_idx - AOI_FilterSignal/Logic - #49, #88, #96, *#9, *#96			
_n0	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_n0 - AOI_FilterSignal/Logic - *#26			
_n1	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_n1 - AOI_FilterSignal/Logic - *#27			
_n2	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_n2 - AOI_FilterSignal/Logic - #75, #77, #77, *#29			
_nBufferUsage	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_nBufferUsage - AOI_FilterSignal/Logic - #50, #50, #51, *#10, *#50			
_nM	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_nM - AOI_FilterSignal/Logic - #16			
_nM - AOI_FilterSignal/Prescan - #4, #6, #8, *#2, *#4, *#6, *#8			
_nN	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
_nN - AOI_FilterSignal/Logic - #29, #30, #32, #34, #36, #39, #50, #51, #72, #86, #88, #96, *#28			

<b>_nTmp</b>	0	DINT	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_nTmp - AOI_FilterSignal/Logic - #16, #28, *#13</i>			
<i>_nTmp - AOI_FilterSignal/Prescan - #4, #4, #6, #8, #8, *#3, *#5, *#7</i>			
<b>_rBufSorted</b>		REAL[31]	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_rBufSorted - AOI_FilterSignal/Logic - #75, #77, #77, *#73</i>			
<i>_rBufSorted - AOI_FilterSignal/Prescan - #5</i>			
<i>_rBufSorted[0] - AOI_FilterSignal/Logic - *#72</i>			
<b>_rRawBuf</b>		REAL[31]	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_rRawBuf - AOI_FilterSignal/Logic - #88, *#49</i>			
<i>_rRawBuf - AOI_FilterSignal/Prescan - #3</i>			
<i>_rRawBuf[0] - AOI_FilterSignal/Logic - #72</i>			
<b>_rSum</b>	0.0	REAL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_rSum - AOI_FilterSignal/Logic - #88, #90, *#85, *#88</i>			
<b>_rSumW</b>	0.0	REAL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_rSumW - AOI_FilterSignal/Logic - #40, #90, *#36, *#38, *#40</i>			
<b>_rWs</b>		REAL[31]	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_rWs - AOI_FilterSignal/Logic - #40, #87</i>			
<i>_rWs - AOI_FilterSignal/Prescan - #2</i>			
<i>_rWs[0] - AOI_FilterSignal/Logic - *#32</i>			
<b>_srt</b>		CONTROL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_srt - AOI_FilterSignal/Logic - *#73</i>			
<b>_srt.LEN</b>	0	DINT	
<i>_srt.LEN - AOI_FilterSignal/Logic - *#34</i>			
<b>_srt.POS</b>	0	DINT	
<i>_srt.POS - AOI_FilterSignal/Logic - *#33</i>			
<b>_w</b>	0.0	REAL	AOI_FilterSignal
Usage:	Local Tag		
External Access:	None		
<i>_w - AOI_FilterSignal/Logic - #88, *#35, *#87</i>			

**AOI\_FilterSignal Instruction Definition - Logic Routine**

test:Add-On Instructions:AOI\_FilterSignal:Logic

Total number of lines in routine: 97

Data Context: AOI\_FilterSignal &lt;definition&gt;

```

if ER and not ( _bFirstScan or In_bInitialize ) then
    tnd();
end_if;

// Initialization
if In_bInitialize or _bFirstScan then
    ER := 0;
    Err := 0;
    _idx := 0;
    _nBufferUsage := 0;

// Validate parameters
_nTmp := Par_nWindowLeft + Par_nWindowsRight + 1;
if Par_nWindowLeft < 0 then Err := 1; end_if;
if Par_nWindowsRight < 0 then Err := 2; end_if;
if _nTmp > _nM then Err := 3; end_if;
if Par_eType <> 0 and Par_eType <> 1 then Err := 4; end_if;
if Err <> 0 then
    ER := 1;
    tnd();
end_if;

// Store parameters ( make filter behaviour independent from outside
// parameters' changes ).
_bW := Par_bUseWeights;
_n0 := Par_nWindowLeft;
_n1 := Par_nWindowsRight;
_nN := _nTmp;
_n2 := _nN / 2;
_bOdd := ( _nN mod 2 ) = 1; // neighbourhood window composed of a odd number of points
_eType := Par_eType;
cop( Par_rWeights[0], _rWs[0], _nN );
_srt.POS := 0;
_srt.LEN := _nN;
_w := 1.0;
_rSumW := _nN;
if _bW then
    _rSumW := 0.0;
    for _i := 0 to _nN - 1 do
        _rSumW := _rSumW + _rWs[_i];
    end_for;
end_if;

end_if;
_bFirstScan := 0;
In_bInitialize := 0;

// While buffering data the filter isn't applied
_rRawBuff[ _idx ] := In_rRawSignal;
if _nBufferUsage < _nN then _nBufferUsage := _nBufferUsage + 1; end_if;
Out_Buffering := _nBufferUsage <> _nN;

if not Out_Buffering then
    case _eType of
        0: // MEDIAN
        // Very simple implementation but not very efficient.

```

**AOI\_FilterSignal Instruction Definition - Logic Routine**

test:Add-On Instructions:AOI\_FilterSignal:Logic

Total number of lines in routine: 97

Data Context: AOI\_FilterSignal &lt;definition&gt;

```

// It's ok here because a single filtered value is
// calculated at every scan of the plc.
// I'm using firmware instruction SRT to
// sort the buffer array, because it's faster
// than anything that could be implemented in
// structure text. In case SRT isn't available in
// other PLCs, a sorting algorithm could be implemented.
// For example, in file sort/sort.c of the Gnu Scientific
// Library (ftp://ftp.gnu.org/gnu/gsl/) there is an
// implementation of the heap sort algorithm.
// TODO: replace the SRT instruction with such GLS
// implementation so that this code can be reused
// in any PLC.

cop( _rRawBuf[0], _rBufSorted[0], _nN );
SRT( _rBufSorted, 0, _srt );
if _bOdd then
    Out_rFilteredSignalDelayed := _rBufSorted[_n2];
else
    Out_rFilteredSignalDelayed := 0.5 * ( _rBufSorted[_n2-1] + _rBufSorted[_n2] );
end_if;

1: // MEAN
// Note: instruction MAVE could be used here, but I find
//       the following code clearer. Moreover the code
//       can be used in any PLC which doesn't include the
//       MAVE instruction.

_rSum := 0.0;
for _i := 0 to _nN - 1 do
    if _bW then _w := _rWs[_i]; end_if;
    _rSum := _rSum + _w * _rRawBuf[ ( _idx + 1 + _i ) mod _nN ];
end_for;
Out_rFilteredSignalDelayed := _rSum / _rSumW;

end_case;

end_if;

_idx := ( _idx + 1 ) mod _nN;

```

**AOI\_FilterSignal Instruction Definition - Prescan Routine**

test:Add-On Instructions:AOI\_FilterSignal:Prescan

Total number of lines in routine: 15

Data Context: AOI\_FilterSignal &lt;definition&gt;


```
// Calculate _nN = minimum between size( Par_rWeights ), size(_rWs), size(_rBuf) and size(_rBufSorted)
size( _rWs, 0, _nM );
size( _rRawBuf, 0, _nTmp );
if _nTmp < _nM then _nM := _nTmp; end_if;
size( _rBufSorted, 0, _nTmp );
if _nTmp < _nM then _nM := _nTmp; end_if;
size( Par_rWeights, 0, _nTmp );
if _nTmp < _nM then _nM := _nTmp; end_if;

// Used for initialization
_bFirstScan := 1;


// Data must be buffered at the beginning
Out_Buffering := 1;
```

**Available Languages**

Relay Ladder



Function Block



Structured Text

```
AOI_Floor(In_X, Out_Ret);
```

**Parameters**

Required	Name	Data Type	Usage	Description
X	AOI_Floor	AOI_Floor	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_X	REAL	Input	
X	Out_Ret	REAL	InOut	

**Extended Description****Execution**

Condition	Description
EnableIn is true	

**Revision v1.1 Notes**

**AOI\_Floor Instruction Definition - Parameter Listing**

test:Add-On Instructions:AOI\_Floor

Data Type Size: 16 byte (s)

Data Context: AOI\_Floor &lt;definition&gt;

Page 129



07-04-2020 15:25:29

LogixProgram\_L27\_Ver4.ACD

Name	Default	Data Type	Scope
<b>In_X</b>	0.0	REAL	AOI_Floor
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_X - AOI_Floor/Logic - #1, #2</i>			
<b>Out_Ret</b>	??	REAL	AOI_Floor
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Out_Ret - AOI_Floor/Logic - *#7</i>			

Name	Default	Data Type	Scope
_n	0	DINT	AOI_Floor
Usage:	Local Tag		
External Access:	None		
_n - AOI_Floor/Logic - #2 #3, #5, *#1			
_ret	0.0	REAL	AOI_Floor
Usage:	Local Tag		
External Access:	None		
_ret - AOI_Floor/Logic - #7, *#3, *#5			

```
_n := In_X;  
if In_X < _n then  
    _ret := _n - 1;  
else  
    _ret := _n;  
end_if;  
Out_Ret := _ret;
```

**Available Languages****Relay Ladder****Function Block****Structured Text**

```
AOI_Max(In_Array, In_MainIndex, In_N0, In_N1);
```

**Parameters**

Required	Name	Data Type	Usage	Description
X	AOI_Max	AOI_Max	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_Array	REAL	InOut	
X	In_MainIndex	DINT	Input	
X	In_N0	INT	Input	
X	In_N1	INT	Input	
	Out_Max	REAL	Output	
	Err	DINT	Output	
	ER	BOOL	Output	
	_i	DINT	Input	

**Extended Description****Execution****Condition      Description**

EnableIn is true

Prescan

**Revision v1.0 Notes**

Name	Default	Data Type	Scope
<b>ER</b>	0	BOOL	AOI_Max
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>ER - AOI_Max/Logic - *#1, *#10</i>			
<b>Err</b>	0	DINT	AOI_Max
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Err - AOI_Max/Logic - #9, *#2, *#5, *#6, *#7, *#8</i>			
<b>In_Array</b>		REAL[1000,2]	AOI_Max
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>In_Array - AOI_Max/Logic - #15, #17</i>			
<i>In_Array - AOI_Max/Prescan - #1, #2</i>			
<b>In_MainIndex</b>	0	DINT	AOI_Max
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_MainIndex - AOI_Max/Logic - #15, #17, #5, #5</i>			
<b>In_N0</b>	0	INT	AOI_Max
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_N0 - AOI_Max/Logic - #15, #16, #6, #6, #8</i>			
<b>In_N1</b>	0	INT	AOI_Max
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_N1 - AOI_Max/Logic - #16, #7, #7, #8</i>			
<b>Out_Max</b>	0.0	REAL	AOI_Max
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Out_Max - AOI_Max/Logic - *#20</i>			
<b>_i</b>	0	DINT	AOI_Max
Usage:	Input Parameter		
Required:	No		
Visible:	No		
External Access:	Read/Write		
<i>_i - AOI_Max/Logic - #17, *#16</i>			

Name	Default	Data Type	Scope
_max	0.0	REAL	AOI_Max
Usage:	Local Tag		
External Access:	None		
_max - AOI_Max/Logic - #18, #20, *#15, *#18			
_maxDim0	0	DINT	AOI_Max
Usage:	Local Tag		
External Access:	None		
_maxDim0 - AOI_Max/Logic - #6, #7			
_maxDim0 - AOI_Max/Prescan - *#1			
_maxDim1	0	DINT	AOI_Max
Usage:	Local Tag		
External Access:	Read/Write		
_maxDim1 - AOI_Max/Logic - #5			
_maxDim1 - AOI_Max/Prescan - *#2			
_tmp	0.0	REAL	AOI_Max
Usage:	Local Tag		
External Access:	None		
_tmp - AOI_Max/Logic - #18, #18, *#17			

```
ER := 0;
Err := 0;

// Some checks
if In_MainIndex < 0 or _maxDim1 <= In_MainIndex then Err := 1; end_if;
if In_N0 < 0          or _maxDim0 <= In_N0      then Err := 2; end_if;
if In_N1 < 0          or _maxDim0 <= In_N1      then Err := 3; end_if;
if In_N0 > In_N1                                then Err := 4; end_if;
if Err <> 0 then
    ER := 1;
    tnd();
end_if;

// Find maximum
_max := In_Array[ In_N0, In_MainIndex ];
for _i := In_N0 + 1 to In_N1 do
    _tmp := In_Array[ _i, In_MainIndex ];
    if _max < _tmp then _max := _tmp; end_if;
end_for;
Out_Max := _max;
```


```
size( In_Array, 0, _maxDim0 ); // Dimension 0 is the most right index
size( In_Array, 1, _maxDim1 );
```

## Available Languages

### Relay Ladder

AOI_Mean2D	
AOI_Mean2D	?
In_Array	?
In_MainIndex	?
	??
In_N0	?
	??
In_N1	?
	??
Out_Mean	?
Err	?

### Function Block



### Structured Text

```
AOI_Mean2D(In_Array, In_MainIndex, In_N0, In_N1);
```

## Parameters

Required	Name	Data Type	Usage	Description
X	AOI_Mean2D	AOI_Mean2D	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_Array	REAL	InOut	
X	In_MainIndex	DINT	Input	
X	In_N0	INT	Input	
X	In_N1	INT	Input	
	Out_Mean	REAL	Output	
	Err	DINT	Output	
	ER	BOOL	Output	
	_i	DINT	Input	

## Extended Description

## Execution

### Condition      Description

EnableIn is true

Prescan

## Revision v1.0 Notes

Name	Default	Data Type	Scope
<b>ER</b>	0	BOOL	AOI_Mean2D
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>ER - AOI_Mean2D/Logic - *#1, *#10</i>			
<b>Err</b>	0	DINT	AOI_Mean2D
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Err - AOI_Mean2D/Logic - #9, *#2, *#5, *#6, *#7, *#8</i>			
<b>In_Array</b>		REAL[1000,2]	AOI_Mean2D
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>In_Array - AOI_Mean2D/Logic - #17</i>			
<i>In_Array - AOI_Mean2D/Prescan - #1, #2</i>			
<b>In_MainIndex</b>	0	DINT	AOI_Mean2D
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_MainIndex - AOI_Mean2D/Logic - #17, #5, #5</i>			
<b>In_N0</b>	0	INT	AOI_Mean2D
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_N0 - AOI_Mean2D/Logic - #16, #19, #6, #6, #8</i>			
<b>In_N1</b>	0	INT	AOI_Mean2D
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_N1 - AOI_Mean2D/Logic - #16, #19, #7, #7, #8</i>			
<b>Out_Mean</b>	0.0	REAL	AOI_Mean2D
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Out_Mean - AOI_Mean2D/Logic - *#19</i>			
<b>_i</b>	0	DINT	AOI_Mean2D
Usage:	Input Parameter		
Required:	No		
Visible:	No		
External Access:	Read/Write		
<i>_i - AOI_Mean2D/Logic - #17, *#16</i>			

Name	Default	Data Type	Scope
_maxDim0	0	DINT	AOI_Mean2D
Usage:	Local Tag		
External Access:	None		
<i>_maxDim0 - AOI_Mean2D/Logic - #6, #7</i>			
<i>_maxDim0 - AOI_Mean2D/Prescan - *#1</i>			
_maxDim1	0	DINT	AOI_Mean2D
Usage:	Local Tag		
External Access:	Read/Write		
<i>_maxDim1 - AOI_Mean2D/Logic - #5</i>			
<i>_maxDim1 - AOI_Mean2D/Prescan - *#2</i>			
_sum	0.0	REAL	AOI_Mean2D
Usage:	Local Tag		
External Access:	None		
<i>_sum - AOI_Mean2D/Logic - #17, #19, *#15, *#17</i>			

```
ER := 0;
Err := 0;

// Some checks
if In_MainIndex < 0 or _maxDim1 <= In_MainIndex then Err := 1; end_if;
if In_N0 < 0      or _maxDim0 <= In_N0      then Err := 2; end_if;
if In_N1 < 0      or _maxDim0 <= In_N1      then Err := 3; end_if;
if In_N0 > In_N1                         then Err := 4; end_if;
if Err <> 0 then
    ER := 1;
    tnd();
end_if;

// Find mean
_sum := 0.0;
for _i := In_N0 to In_N1 do
    _sum := _sum + In_Array[ _i, In_MainIndex ];
end_for;
Out_Mean := _sum / ( In_N1 - In_N0 + 1 );
```


```
size( In_Array, 0, _maxDim0 ); // Dimension 0 is the most right index
size( In_Array, 1, _maxDim1 );
```

## Available Languages

### Relay Ladder

AOI_Mean50	
AOI_Mean50	?
In_Array	?
In_N0	?
	??
In_N1	?
	??
Out_Mean	?
Err	?

### Function Block



### Structured Text

```
AOI_Mean50(In_Array, In_N0, In_N1);
```

## Parameters

Required	Name	Data Type	Usage	Description
X	AOI_Mean50	AOI_Mean50	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_Array	REAL	InOut	
X	In_N0	INT	Input	
X	In_N1	INT	Input	
	Out_Mean	REAL	Output	
	Err	DINT	Output	
	ER	BOOL	Output	
	_i	DINT	Input	

## Extended Description

## Execution

Condition	Description
-----------	-------------

EnableIn is true  
Prescan

## Revision v1.0 Notes

Name	Default	Data Type	Scope
<b>ER</b>	0	BOOL	AOI_Mean50
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>ER - AOI_Mean50/Logic - *#1, *#9</i>			
<b>Err</b>	0	DINT	AOI_Mean50
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Err - AOI_Mean50/Logic - #8, *#2, *#5, *#6, *#7</i>			
<b>In_Array</b>		REAL[50]	AOI_Mean50
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>In_Array - AOI_Mean50/Logic - #16</i>			
<i>In_Array - AOI_Mean50/Prescan - #1</i>			
<b>In_N0</b>	0	INT	AOI_Mean50
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_N0 - AOI_Mean50/Logic - #15, #18, #5, #5, #7</i>			
<b>In_N1</b>	0	INT	AOI_Mean50
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_N1 - AOI_Mean50/Logic - #15, #18, #6, #6, #7</i>			
<b>Out_Mean</b>	0.0	REAL	AOI_Mean50
Usage:	Output Parameter		
Required:	No		
Visible:	Yes		
External Access:	Read/Write		
<i>Out_Mean - AOI_Mean50/Logic - *#18</i>			
<b>_i</b>	0	DINT	AOI_Mean50
Usage:	Input Parameter		
Required:	No		
Visible:	No		
External Access:	Read/Write		
<i>_i - AOI_Mean50/Logic - #16, *#15</i>			

Name	Default	Data Type	Scope
_maxDim0	0	DINT	AOI_Mean50
Usage: External Access:	Local Tag None		
	_maxDim0 - AOI_Mean50/Logic - #5, #6		
	_maxDim0 - AOI_Mean50/Prescan - *#1		
_sum	0.0	REAL	AOI_Mean50
Usage: External Access:	Local Tag None		
	_sum - AOI_Mean50/Logic - #16, #18, *#14, *#16		

```
ER := 0;
Err := 0;

// Some checks
if In_N0 < 0      or _maxDim0 <= In_N0      then Err := 2; end_if;
if In_N1 < 0      or _maxDim0 <= In_N1      then Err := 3; end_if;
if In_N0 > In_N1      then Err := 4; end_if;
if Err <> 0 then
    ER := 1;
    tnd();
end_if;

// Find mean
_sum := 0.0;
for _i := In_N0 to In_N1 do
    _sum := _sum + In_Array[ _i];
end_for;
Out_Mean := _sum / ( In_N1 - In_N0 + 1 );
```

---


```
size( In_Array, 0, _maxDim0 ); // Dimension 0 is the most right index
```

**Available Languages**

Relay Ladder



Function Block



Structured Text

```
AOI_Round(In_X, Out_Ret);
```

**Parameters**

Required	Name	Data Type	Usage	Description
X	AOI_Round	AOI_Round	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_X	REAL	Input	
X	Out_Ret	REAL	InOut	

**Extended Description****Execution****Condition      Description**

EnableIn is true

**Revision v1.0 Notes**

Name	Default	Data Type	Scope
<b>In_X</b>	0.0	REAL	AOI_Round
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_X - AOI_Round/Logic - #1</i>			
<b>Out_Ret</b>	??	REAL	AOI_Round
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Out_Ret - AOI_Round/Logic - *#2</i>			

Name	Default	Data Type	Scope
_n	0	DINT	AOI_Round


Usage: Local Tag  
External Access: None  
*\_n - AOI\_Round/Logic - #2, \*#I*

---


```
_n := In_X;  
Out_Ret := _n;
```

## Available Languages

Relay Ladder



Function Block



Structured Text

```
AOI_Trunc(In_X, Out_Ret);
```

## Parameters

Required	Name	Data Type	Usage	Description
X	AOI_Trunc	AOI_Trunc	InOut	
	EnableIn	BOOL	Input	
	EnableOut	BOOL	Output	
X	In_X	REAL	Input	
X	Out_Ret	REAL	InOut	

## Extended Description

## Execution

Condition	Description
EnableIn is true	

## Revision v1.0 Notes

Name	Default	Data Type	Scope
<b>In_X</b>	0.0	REAL	AOI_Trunc
Usage:	Input Parameter		
Required:	Yes		
Visible:	Yes		
External Access:	Read/Write		
<i>In_X - AOI_Trunc/Logic - #7</i>			
<b>Out_Ret</b>	??	REAL	AOI_Trunc
Usage:	InOut Parameter		
Required:	Yes		
Visible:	Yes		
Constant	No		
<i>Out_Ret - AOI_Trunc/Logic - *#7</i>			

Name	Default	Data Type	Scope
<b>No Tags Exist</b>			

```
// TODO: Logix has a trunc function but other PLCs may not.  
// Replace the native function in Logix with an  
// implementation which will work with other PLC.  
// Note that AOI_Trunc and AOI_Ceil (which depends  
// on AOI_Trunc ) may not be needed in the project.
```

```
Out_Ret := trunc( In_X );
```

Data type Name: STRING

Description:

Size: 88 byte(s)

Name	Value	Data Type	Style
<b>LEN</b> External Access:	Read/Write	DINT	Decimal
<b>DATA</b> External Access:	Read/Write	SINT[82]	ASCII

Data type Name: AOI\_ARRAYSHIFT

Description:

Size: 4 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		

LogixProgram\_L27\_Ver4.ACD

Data type Name: AOI\_Ceil

Description:

Size: 24 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_X</b>		REAL	Float
External Access:	Read/Write		

Data type Name: AOI\_FillArrayInt

Description:

Size: 12 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>Source</b>		DINT	Decimal
External Access:	Read/Write		
<b>Len</b>		DINT	Decimal
External Access:	Read/Write		

Data type Name: AOI\_FillArrayReal

Description:

Size: 12 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>Source</b>		DINT	Decimal
External Access:	Read/Write		
<b>Len</b>		DINT	Decimal
External Access:	Read/Write		

Data type Name: AOI\_FillArrayReal50

Description:

Size: 12 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>Source</b>		DINT	Decimal
External Access:	Read/Write		
<b>Len</b>		DINT	Decimal
External Access:	Read/Write		

Data type Name: AOI\_FillArray\_Real2D

Description:

Size: 12 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>Source</b>		DINT	Decimal
External Access:	Read/Write		
<b>Len</b>		DINT	Decimal
External Access:	Read/Write		

LogixProgram\_L27\_Ver4.ACD

Data type Name: AOI\_FilterSignal

Description:

Size: 468 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_rRawSignal</b>		REAL	Float
External Access:	Read/Write		
<b>In_bInitialize</b>		BOOL	Decimal
External Access:	Read/Write		
<b>Par_eType</b>		DINT	Decimal
0: MEDIAN; 1: MEAN			
External Access:	Read/Write		
<b>Par_nWindowLeft</b>		DINT	Decimal
External Access:	Read/Write		
<b>Par_nWindowsRight</b>		DINT	Decimal
External Access:	Read/Write		
<b>Par_bUseWeights</b>		BOOL	Decimal
Used for mean filtering			
External Access:	Read/Write		
<b>Out_Buffering</b>		BOOL	Decimal
External Access:	Read/Write		
<b>Out_rFilteredSignalDelayed</b>		REAL	Float
External Access:	Read Only		
<b>ER</b>		BOOL	Decimal
External Access:	Read Only		
<b>Err</b>		DINT	Decimal
External Access:	Read/Write		

LogixProgram\_L27\_Ver4.ACD

Data type Name: AOI\_Floor

Description:

Size: 16 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_X</b>		REAL	Float
External Access:	Read/Write		

LogixProgram\_L27\_Ver4.ACD

Data type Name: AOI\_Max

Description:

Size: 40 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_MainIndex</b>		DINT	Decimal
External Access:	Read/Write		
<b>In_N0</b>		INT	Decimal
External Access:	Read/Write		
<b>In_N1</b>		INT	Decimal
External Access:	Read/Write		
<b>Out_Max</b>		REAL	Float
External Access:	Read/Write		
<b>Err</b>		DINT	Decimal
External Access:	Read/Write		
<b>ER</b>		BOOL	Decimal
External Access:	Read/Write		
<b>i</b>		DINT	Decimal
External Access:	Read/Write		

Data type Name: AOI\_Mean2D

Description:

Size: 36 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_MainIndex</b>		DINT	Decimal
External Access:	Read/Write		
<b>In_N0</b>		INT	Decimal
External Access:	Read/Write		
<b>In_N1</b>		INT	Decimal
External Access:	Read/Write		
<b>Out_Mean</b>		REAL	Float
External Access:	Read/Write		
<b>Err</b>		DINT	Decimal
External Access:	Read/Write		
<b>ER</b>		BOOL	Decimal
External Access:	Read/Write		
<b>i</b>		DINT	Decimal
External Access:	Read/Write		

Data type Name: AOI\_Mean50

Description:

Size: 32 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_N0</b>		INT	Decimal
External Access:	Read/Write		
<b>In_N1</b>		INT	Decimal
External Access:	Read/Write		
<b>Out_Mean</b>		REAL	Float
External Access:	Read/Write		
<b>Err</b>		DINT	Decimal
External Access:	Read/Write		
<b>ER</b>		BOOL	Decimal
External Access:	Read/Write		
<b>i</b>		DINT	Decimal
External Access:	Read/Write		

LogixProgram\_L27\_Ver4.ACD

Data type Name: AOI\_Round

Description:

Size: 12 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_X</b>		REAL	Float
External Access:	Read/Write		

Data type Name: AOI\_Trunc

Description:

Size: 8 byte(s)

Name	Value	Data Type	Style
<b>EnableIn</b>		BOOL	Decimal
Enable Input - System Defined Parameter			
External Access:	Read Only		
<b>EnableOut</b>		BOOL	Decimal
Enable Output - System Defined Parameter			
External Access:	Read Only		
<b>In_X</b>		REAL	Float
External Access:	Read/Write		

**1769 Bus : Local Modules****Local: [0] 1769-L27ERM-QBFC1B test**

Type:	1769-L27ERM-QBFC1B CompactLogix™ 5370 Controller	Parent:	Local
Vendor:	Rockwell Automation/Allen-Bradley	Vendor ID:	1
Slot:	0	Electronic Keying:	Disabled
Revision:	32.13	Status:	Standby
Module Fault:	Offline	Inhibit Flag	Off

**Embedded I/O : Local Modules****Local: [1] Embedded Discrete\_IO**

Type:	Embedded 16 Point 24V DC Sink/Source Input /16 Point 24V DC Source Output	Parent:	Local
Vendor:	Rockwell Automation/Allen-Bradley	Vendor ID:	1
Slot:	1	Electronic Keying:	Compatible Keying
Revision:	3.1	Status:	Standby
Module Fault:	Offline	Inhibit Flag	Off
Use Unicast:	No		

**Module Defined Configuration Tag**

Value	Data Type
Local:1:C	AB:Embedded_DiscreteIO1:C:0
.Filter0OffOn_0	BOOL
.Filter0OffOn_1	BOOL
.Filter0OffOn_2	BOOL
.Filter0OffOn_3	BOOL
.Filter0OnOff_4	BOOL
.Filter0OnOff_5	BOOL
.Filter0OnOff_6	BOOL
.Filter0OnOff_7	BOOL
.Filter1OffOn_0	BOOL
.Filter1OffOn_1	BOOL
.Filter1OffOn_2	BOOL
.Filter1OffOn_3	BOOL
.Filter1OnOff_4	BOOL
.Filter1OnOff_5	BOOL
.Filter1OnOff_6	BOOL
.Filter1OnOff_7	BOOL

**Local: [2] Embedded Analog\_IO**

Type:	Embedded 4 Channel Universal Analog Input/2 Channel Analog Output	Parent:	Local
Vendor:	Rockwell Automation/Allen-Bradley	Vendor ID:	1
Slot:	2	Electronic Keying:	Compatible Keying
Revision:	1.1	Status:	Standby
Module Fault:	Offline	Inhibit Flag	Off
Use Unicast:	No		

**Module Defined Configuration Tag**

Value	Data Type
Local:2:C	AB:Embedded_AnalogIO1:C:0
.RealTimeSample	INT
.TimestampEn	BOOL
.Ch0InputFilter_0	BOOL
.Ch0InputFilter_1	BOOL
.Ch0InputFilter_2	BOOL
.Ch0InputFilter_3	BOOL

LogixProgram\_L27\_Ver4.ACD

.Ch0DisableCJComp	1	BOOL
ensation		
.Ch0OpenWire_5	0	BOOL
.Ch0OpenWire_6	0	BOOL
.Ch0OpenWireEn	0	BOOL
.Ch0AlarmLatchEn	0	BOOL
.Ch0AlarmEn	1	BOOL
.Ch0InputEn	1	BOOL
.Ch0InputRange	4	SINT
	1	BOOL
.Ch0InputDataFormat		
_0	0	BOOL
.Ch0InputDataFormat		
_1	0	BOOL
.Ch0InputDataFormat		
_2	0	BOOL
.Ch0TempUnits	0	BOOL
	0	BOOL
.Ch0RTD_Resistan		
ceMeasType_6	0	BOOL
.Ch0RTD_Resistan		
ceMeasType_7		
.Ch0HAlarmLimit	2000	INT
.Ch0LAlarmLimit	500	INT
	0	INT
.Ch0AlarmDeadband		
.Ch1InputFilter_0	1	BOOL
.Ch1InputFilter_1	0	BOOL
.Ch1InputFilter_2	1	BOOL
.Ch1InputFilter_3	0	BOOL
	1	BOOL
.Ch1DisableCJComp		
ensation		
.Ch1OpenWire_5	0	BOOL
.Ch1OpenWire_6	0	BOOL
.Ch1OpenWireEn	0	BOOL
.Ch1AlarmLatchEn	0	BOOL
.Ch1AlarmEn	1	BOOL
.Ch1InputEn	1	BOOL
.Ch1InputRange	4	SINT
	1	BOOL
.Ch1InputDataFormat		
_0	0	BOOL
.Ch1InputDataFormat		
_1	0	BOOL
.Ch1InputDataFormat		
_2	0	BOOL
.Ch1TempUnits	0	BOOL
.Ch1HAlarmLimit	5250	INT
.Ch1LAlarmLimit	500	INT
	0	INT
.Ch1AlarmDeadband		
.Ch2InputFilter_0	1	BOOL
.Ch2InputFilter_1	0	BOOL
.Ch2InputFilter_2	1	BOOL
.Ch2InputFilter_3	0	BOOL

LogixProgram\_L27\_Ver4.ACD

.Ch2DisableCJComp	1	BOOL
ensation		
.Ch2OpenWire_5	0	BOOL
.Ch2OpenWire_6	0	BOOL
.Ch2OpenWireEn	0	BOOL
.Ch2AlarmLatchEn	1	BOOL
.Ch2AlarmEn	1	BOOL
.Ch2InputEn	1	BOOL
.Ch2InputRange	3	SINT
	1	BOOL
.Ch2InputDataFormat		
_0	0	BOOL
.Ch2InputDataFormat		
_1	0	BOOL
.Ch2InputDataFormat		
_2	0	BOOL
.Ch2TempUnits	0	BOOL
	0	BOOL
.Ch2RTD_Resistan		
ceMeasType_6	0	BOOL
.Ch2RTD_Resistan		
ceMeasType_7		
.Ch2HAlarmLimit	5600	INT
.Ch2LAlarmLimit	3200	INT
	0	INT
.Ch2AlarmDeadband		
.Ch3InputFilter_0	1	BOOL
.Ch3InputFilter_1	0	BOOL
.Ch3InputFilter_2	1	BOOL
.Ch3InputFilter_3	0	BOOL
	1	BOOL
.Ch3DisableCJComp		
ensation		
.Ch3OpenWire_5	0	BOOL
.Ch3OpenWire_6	0	BOOL
.Ch3OpenWireEn	0	BOOL
.Ch3AlarmLatchEn	0	BOOL
.Ch3AlarmEn	0	BOOL
.Ch3InputEn	1	BOOL
.Ch3InputRange	1	SINT
	1	BOOL
.Ch3InputDataFormat		
_0	0	BOOL
.Ch3InputDataFormat		
_1	0	BOOL
.Ch3InputDataFormat		
_2		
.Ch3TempUnits	0	BOOL
.Ch3HAlarmLimit	2	INT
.Ch3LAlarmLimit	0	INT
	0	INT
.Ch3AlarmDeadband		
.CJTempUnits	0	BOOL
.CJWeightedProfile	0	BOOL
	0	BOOL
.CyclicCalibrationEn		

LogixProgram\_L27\_Ver4.ACD

.UpdateCJCompensationEn	0	BOOL
.Ch0LimitAlarmLatchedEn	0	BOOL
.Ch0OutputEn	1	BOOL
	3	SINT
.Ch0OutputRangeType	1	SINT
.Ch0OutputDataFormat		
.Ch0HClampValue	21000	INT
.Ch0LClampValue	3200	INT
	0	BOOL
.Ch1LimitAlarmLatchedEn		
.Ch1OutputEn	1	BOOL
	3	SINT
.Ch1OutputRangeType	1	SINT
.Ch1OutputDataFormat		
.Ch1HClampValue	21000	INT
.Ch1LClampValue	3200	INT

**Local: [3] Embedded Counters**

Type:	Embedded 4 Channel High Speed Counter	Parent:	Local
Vendor:	Rockwell Automation/Allen-Bradley	Vendor ID:	1
Slot:	3	Electronic Keying:	Compatible Keying
Revision:	1.1	Status:	Standby
Module Fault:	Offline	Inhibit Flag	Off
Use Unicast:	No		

Module Defined Configuration Tag	Value	Data Type
Local:3:C		AB:Embedded_HSC1:C:0
.Config0	2#0000_0001_0000_0000	INT
	0	BOOL
.OverCurrentLatchOff		
.CtrReset	0	BOOL
.ProgToFaultEn	0	BOOL
	1	BOOL
.NumberOfCounters_0	0	BOOL
.NumberOfCounters_1	0	BOOL
.Ctr0ResetEn	0	BOOL
.Ctr1ResetEn	0	BOOL
.Ctr2ResetEn	0	BOOL
.Ctr3ResetEn	0	BOOL
.Filter	2#0000_0000_0000_0000	INT
.FilterA0_0	0	BOOL
.FilterA0_1	0	BOOL
.FilterB0_0	0	BOOL
.FilterB0_1	0	BOOL
.FilterZ0_0	0	BOOL
.FilterZ0_1	0	BOOL
.FilterA1_0	0	BOOL

LogixProgram\_L27\_Ver4.ACD

.FilterA1_1	0	BOOL
.FilterB1_0	0	BOOL
.FilterB1_1	0	BOOL
.FilterZ1_0	0	BOOL
.FilterZ1_1	0	BOOL
.Ctr0MaxCount	2147483647	DINT
.Ctr0MinCount	-2147483648	DINT
.Ctr0Preset	0	DINT
.Ctr0Hysteresis	0	INT
.Ctr0Scaler	1	INT
	10	INT
.Ctr0CyclicRateUpdateTime		
.Ctr0Config	2#0000_0000	SINT
	0	BOOL
.Ctr0OperationMode_0	0	BOOL
.Ctr0OperationMode_1	0	BOOL
.Ctr0OperationMode_2	0	BOOL
.Ctr0StoreOnRisingZ		
.Ctr0HoldOnZ	0	BOOL
	0	BOOL
.Ctr0PresetOnRisingZ		
.Ctr0Linear	0	BOOL
.Ctr1MaxCount	2147483647	DINT
.Ctr1MinCount	-2147483648	DINT
.Ctr1Preset	0	DINT
.Ctr1Hysteresis	0	INT
.Ctr1Scaler	1	INT
	10	INT
.Ctr1CyclicRateUpdateTime		
.Ctr1Config	2#0000_0000	SINT
	0	BOOL
.Ctr1OperationMode_0	0	BOOL
.Ctr1OperationMode_1	0	BOOL
.Ctr1OperationMode_2	0	BOOL
.Ctr1StoreOnRisingZ		
.Ctr1HoldOnZ	0	BOOL
	0	BOOL
.Ctr1PresetOnRisingZ		
.Ctr1Linear	0	BOOL
.Ctr2MaxCount	0	DINT
.Ctr2MinCount	0	DINT
.Ctr2Preset	0	DINT
.Ctr2Hysteresis	0	INT
.Ctr2Scaler	0	INT
	0	INT
.Ctr2CyclicRateUpdateTime		
.Ctr2Linear	0	BOOL
.Ctr3MaxCount	0	DINT
.Ctr3MinCount	0	DINT

LogixProgram\_L27\_Ver4.ACD

.Ctr3Preset	0	DINT
.Ctr3Hysteresis	0	INT
.Ctr3Scaler	0	INT
	0	INT
.Ctr3CyclicRateUpdateTime		
.Ctr3Linear	0	BOOL
.Range0To11		AB:1769_HSC1_Range:C:0[12]

LogixProgram\_L27\_Ver4.ACD

<b>test</b>	
Controller Organizer Listing.....	1
Controller Properties Listing.....	3
Tag Listing .....	4
<b>P200ms_p10</b>	
Task Properties Listing.....	30
<b>GUI</b>	
Program Properties Listing.....	31
<b>Handle_GUI</b>	
Structured Text.....	32
Routine Properties Listing.....	36
<b>ValveProfile</b>	
Program Properties Listing.....	37
<b>Main</b>	
Ladder Diagram .....	38
Routine Properties Listing.....	40
<b>P_Sample_p5</b>	
Task Properties Listing.....	41
<b>measArray</b>	
Program Properties Listing.....	42
<b>FirstScan</b>	
Structured Text.....	43
Routine Properties Listing.....	46
<b>Main</b>	
Structured Text.....	47
Routine Properties Listing.....	48
<b>measArrays</b>	
Structured Text.....	49
Routine Properties Listing.....	52
<b>ISR_ValveCtrl</b>	
Program Properties Listing.....	53
<b>Main</b>	
Structured Text.....	54
Routine Properties Listing.....	58
<b>P_Sample_p8</b>	
Task Properties Listing.....	59
<b>Setup</b>	
Program Properties Listing.....	60
<b>InitGUI</b>	
Structured Text.....	61
Routine Properties Listing.....	62
<b>InitISR_RS232</b>	
Structured Text.....	63
Routine Properties Listing.....	64
<b>InitVariables</b>	
Structured Text.....	65
Routine Properties Listing.....	66
<b>Main</b>	
Structured Text.....	67
Routine Properties Listing.....	68
<b>calcBreath</b>	
Program Properties Listing.....	69
<b>Main</b>	
Structured Text.....	70
Routine Properties Listing.....	72
<b>initVent</b>	
Program Properties Listing.....	73
<b>initISR_ValveCtrl</b>	
Structured Text.....	74
Routine Properties Listing.....	75
<b>Init_ISR_PressWatchID</b>	
Structured Text.....	76
Routine Properties Listing.....	77
<b>Main</b>	

LogixProgram\_L27\_Ver4.ACD

Structured Text.....	78
Routine Properties Listing.....	79
<b>BreathDetect</b>	
Program Properties Listing.....	80
<b>Main</b>	
Structured Text.....	81
Routine Properties Listing.....	82
<b>measCalc</b>	
Program Properties Listing.....	83
<b>Main</b>	
Structured Text.....	84
Routine Properties Listing.....	89
<b>MainFlow</b>	
Program Properties Listing.....	90
<b>Main</b>	
Structured Text.....	91
Routine Properties Listing.....	93
<b>Add-On Instruction Signature Listing</b>	
<b>Add-On Instructions</b>	
<b>AOI_ARRAYSHIFT</b>	
Instruction Definition .....	95
Parameter Listing .....	96
Local Tag Listing .....	97
Logic Routine.....	98
<b>AOI_Ceil</b>	
Instruction Definition .....	99
Parameter Listing .....	100
Local Tag Listing .....	101
Logic Routine.....	102
<b>AOI_FillArrayInt</b>	
Instruction Definition .....	103
Parameter Listing .....	104
Local Tag Listing .....	105
Logic Routine.....	106
<b>AOI_FillArrayReal</b>	
Instruction Definition .....	107
Parameter Listing .....	108
Local Tag Listing .....	109
Logic Routine.....	110
<b>AOI_FillArrayReal50</b>	
Instruction Definition .....	111
Parameter Listing .....	112
Local Tag Listing .....	113
Logic Routine.....	114
<b>AOI_FillArray_Real2D</b>	
Instruction Definition .....	115
Parameter Listing .....	116
Local Tag Listing .....	117
Logic Routine.....	118
<b>AOI_FilterSignal</b>	
Instruction Definition .....	119
Parameter Listing .....	121
Local Tag Listing .....	123
Logic Routine.....	125
Prescan Routine.....	127
<b>AOI_Floor</b>	
Instruction Definition .....	128
Parameter Listing .....	129
Local Tag Listing .....	130
Logic Routine.....	131
<b>AOI_Max</b>	
Instruction Definition .....	132
Parameter Listing .....	133
Local Tag Listing .....	134

LogixProgram\_L27\_Ver4.ACD

Logic Routine.....	135
Prescan Routine.....	136
<b>AOI_Mean2D</b>	
Instruction Definition .....	137
Parameter Listing .....	138
Local Tag Listing .....	139
Logic Routine.....	140
Prescan Routine.....	141
<b>AOI_Mean50</b>	
Instruction Definition .....	142
Parameter Listing .....	143
Local Tag Listing .....	144
Logic Routine.....	145
Prescan Routine.....	146
<b>AOI_Round</b>	
Instruction Definition .....	147
Parameter Listing .....	148
Local Tag Listing .....	149
Logic Routine.....	150
<b>AOI_Trunc</b>	
Instruction Definition .....	151
Parameter Listing .....	152
Local Tag Listing .....	153
Logic Routine.....	154
<b>Data Types</b>	
Strings .....	155
Add-On-Defined Data Type.....	156
<b>Module Properties</b>	
1769 Bus : Local Modules .....	169
Embedded I/O : Local Modules .....	169